


Bill To: (EasyStand Supplier Only)	
Purchase Order #:	
Quote #:	Account #:
Ordered By:	Date:
Supplier:	
Address:	
City, State, Zip:	
Phone:	
Fax:	

P7000ei EasyStand[®] 7000 Magician•ei.....\$1,842

The Magician•ei comes standard with royal blue or ruby red wrap-around seat and back support with built-in seat belt and chest strap, kneepads, foot plates, 4" locking casters, clear acrylic tray with front pad, and gas spring lift.



Upholstery Color: (royal blue will be selected if no color is chosen)

- Royal Blue (standard).....N/C
- Ruby Red.....N/C

Ship To:	
Name:	
Address:	
Address:	
City, State, Zip:	
Phone:	
Tag For:	

Options:


- P80878 Foot Straps (pair).....\$ 36
- P80864 Secure Foot Option (pair).....\$149
- P80873 Tiny Toes Adapters (pair, requires P80878 or P80864).....\$128
- P80848 Independent Kneepads (replaces standard kneepads).....\$180
- P80998 Angle Adjustable Tray for 7000.....\$293
- P5024 Tray Swing-away Sleeve.....\$ 56
- P5070 Ratchet Handles.....\$ 97
- P80885 Shoulder Harness.....\$149
- P81149 Chest Vest-Small.....\$190
- P81292 Royal Blue Form-to-Fit Upholstery (back & seat).....\$324
- P81293 Ruby Red Form-to-Fit Upholstery (back & seat).....\$324
- P80809 Tool Pouch.....\$ 20
- P81356 Royal Blue Head Support.....\$303
- P81357 Ruby Red Head Support.....\$303
- P80869 Stealth[®] Headrest Ball Adapter (.656").....\$149
- P80870 Whitmyer Headrest Ball Adapter (.75").....\$149
- P80871 Whitmyer Headrest Ball Adapter (.625").....\$149
- P80872 Otto Bock[™] Headrest Ball Adapter (1").....\$149

Client Info:		
Height:	Weight:	Diagnosis:
A:	B:	
Facility:		
Therapist:		

Notes:

Magician•ei Specifications

Approximate Height Range	28"-40"
Weight Limit	45 lbs.
A: Seat Depth Range (seat pivot to hip pivot)	6.5"-11"
B: Seat to Footplate Range (from seat pivot)	3"-15"



A: Measure the distance between the back of the buttocks and the back of the knee.

B: Measure the distance between the foot/shoe and the knee bend.