



Quickie Revolution

EI2D, Effective July 2010

ACCOUNT INFORMATION

Account No: _____ Drop Ship/Ship To: _____
 Date: _____ Name: _____
 P.O. #: _____ Address: _____
 Buyer: _____ City: _____ State: _____
 Marked For: _____ Zip Code: _____ Phone Number: _____

The HCPCS codes provided are based on code verification by the PDAC or our interpretation of the code definitions.
 Proper use of HCPCS codes does not ensure coverage or payment. For coverage information, verify the appropriate payor's coverage policy.

QUICKIE REVOLUTION

033WM19 K0005 Quickie Revolution 75°V \$1,999.00

FRAMES

FRAME COLOR

- 033C2 Black NC
- 033C10 Black Opal NC
- 033C38 Blue Opal NC
- 033C36 Candy Blue NC
- 033C35 Candy Red NC
- 033C128 Copper NC
- 033C6 Forest Green NC
- 033C127 Fuchsia NC
- 033C126 Glow NC
- 033C4 Midnight Purple NC
- 033C5 Pearl Pink NC
- 033C103 Rootbeer NC
- 033C114 Sparkle Silver NC
- 033C113 Titanium Color Paint NC
- 033C11 Toxic Green NC
- 033C8 Yellow NC

CUSTOM FRAME COLOR

- 033C112 Desert Camouflage \$295.00
- 033C43 Mossy Oak Camouflage \$295.00
- 033C42 Stars & Stripes \$295.00
- 033C108 Zebra (Choose above base color) \$295.00
*Std with White unless above base paint is chosen
 NA with Black.*

Frame Width

Seat width 1/2" narrower

- 033W5 14" NC
- 033W6 15" NC
- 033W7 16" NC
- 033W8 17" NC
- 033W9 18" NC
- 033W10 19" NC
- 033W11 E2201 20" \$450.00

Sling Depth

- 033SD7 14" NC
- 033SD8 15" NC
- 033SD9 16" NC
- 033SD10 17" NC
- 033SD11 18" NC

Sling Depth Note: The Quickie Revolution has an additional 1" space from seat sling to backrest (compared to other Quickie chairs). In consideration of this a 1" shorter sling may be ordered.

SEATING

Seat Sling

- 033SL1 Standard NC
- 033SL6 Omit Sling -\$15.00

Cushion

- 033CU1 K0669 2" Cushion \$95.00
- 033CU2 Omit Cushion NC

BACKRESTS

Back Type - Folding

- 033BP26 Straight Post NC
- 033BP25 Push Handle \$55.00
Push handles 8° for 16"-20" back size

Back Height

- 033B83 12" - 16" NC
- 033B94 16" - 20" NC

Back Upholstery

- 033BU1 Standard NC
- 033BU8 Adjustable \$175.00
- 033BU14 Omit Back Upholstery -\$15.00

ARMRESTS

Height Adjustable

NA w /8° or 12° camber or folding aluminum side guards; must be removed to fold the backrest down

- 033AA8 E0973 Adult \$260.00
- 033AM1 Desk Length (10")
- 033AM2 Full Length (14")
- 033AA9 E0973 Kids \$260.00
- 033AM1 Desk Length (10")
- 033AM2 Full Length (14")
- 033AA15 E0973 Padded Swing Away \$175.00

FOOTREST

Footrest (Includes Leg Strap)

- 033F52 K0040 Platform Standard (14" to 17-3/4") \$95.00
 K0038
- 033F51 K0040 Platform Long (16-1/4" to 20") \$95.00
 K0038

Footrest Options

- 033LS8 Omit Leg Strap NC

CASTERS

- 033CW1 3" Micro NC
- 5" Low-profile Polyurethane Std
- 033CT6 5" x 1.5" Semi-Pneumatic \$85.00
Only available for use in bottom fork holes
- 6" Polyurethane NC
- 6" x 1.5" Semi-Pneumatic \$85.00
Only available for use in bottom fork holes

Caster Options

- 033CS5 3/4" Longer Fork Stem Bolt NC
- 033CS9 1-1/2" Longer Fork Stem Bolt NC
- 033CS3 K0073 Caster Pin Locks \$195.00

AXLE AND WHEEL (See Charts on Right)

Axle Bracket

<input type="checkbox"/>	033AL23	Standard	NC
<input type="checkbox"/>	033AL25	+3/4" Spacer	NC
<input type="checkbox"/>	033AL24	+1-1/2" Spacer	NC

Camber Tube

<input type="checkbox"/>	033C11	0°	NC
<input type="checkbox"/>	033C13	4°	NC
<input type="checkbox"/>	033C15	8°	NC
<input type="checkbox"/>	033C17	12°	NC

NA w/ height adjustable armrests

Axle

<input type="checkbox"/>	033AX1	Quick Release Axle	NC
<input type="checkbox"/>	033AX11	Quad Quick Release Axle	\$90.00

Rear Anti-Tip

<input type="checkbox"/>	033AT1	E0971 Rear Anti-Tip	\$135.00
<input type="checkbox"/>	033AT3	Omit Rear Anti-Tip	NC

Rear Wheel Type

<input type="checkbox"/>		22" Mag	NC
<input type="checkbox"/>		22" Spoke	\$200.00
<input type="checkbox"/>		24" Mag	NC
<input type="checkbox"/>		24" Spoke	NC
<input type="checkbox"/>		26" Spoke	\$60.00
<input type="checkbox"/>	033RR25	Omit Rear Wheel & Axle (Intended Size for Set Up)	-\$55.00

Quickie Performance

Contact Built 4 Me for Spinergy color options

<input type="checkbox"/>		24" Spinergy 540 (Clinchers Only)	\$695.00
--------------------------	--	-----------------------------------	----------

Rear Wheel Tire

<input type="checkbox"/>	033RT1	E2211 Pneumatic	NC
		E2212 <i>NA w/ performance wheels</i>	
<input type="checkbox"/>	033RT2	E2213 Pneumatic Airless Insert	\$80.00
		E2211 <i>NA w/ 26" or performance wheels</i>	
<input type="checkbox"/>	033RT3	E2218 Full Polyurethane	NC
		<i>NA w/ 26" or performance wheels</i>	
<input type="checkbox"/>	033RT6	E2211 High Pressure Clincher	\$95.00
		E2212 <i>NA w/ 22" or mag wheels</i>	

Handrim

<input type="checkbox"/>	033HR1	Aluminum Anodize	NC
<input type="checkbox"/>	033HR4	Plastic Coated	\$95.00
<input type="checkbox"/>	033HR7	Omit Handrim	NC

Projections

NA w/ Plastic Coated Handrims

<input type="checkbox"/>	033HR8	8 Oblique	\$240.00
		<i>NA w/ 22" wheels.</i>	

Wheel Lock

<input type="checkbox"/>	033WL6	High Mount Scissor	NC
<input type="checkbox"/>	033WL4	Swing-Away	NC
		<i>NA w/ 8° or 12° camber, height adj armrests, plastic or folding aluminum sideguards</i>	
<input type="checkbox"/>	033WL12	Do Not Mount Wheel Lock	NC
<input type="checkbox"/>	033WL10	Omit Wheel Lock	NC

Rear Seat Height	Camber	Wheel Size and Spacer
14 3/4"	8°, 12°	22"
15 1/4"	0°, 4°	22"
	8°, 12°	24" HP
15 1/2"	8°, 12°	22" + 3/4" Spacer
15 3/4"	0°, 4°	24" HP
	8°, 12°	24"
16"	0°, 4°	22" + 3/4" Spacer
	8°, 12°	24"
16 1/4"	8°, 12°	22" + 1 1/2" Spacer
	0°, 4°	24"
16 1/2"	0°, 4°	24" HP + 3/4" Spacer
	8°, 12°	24" + 3/4" Spacer
16 3/4"	0°, 4°	22" + 1 1/2" Spacer
	8°, 12°	24" HP + 1 1/2" Spacer
	8°, 12°	26"
17"	0°, 4°	24" + 3/4" Spacer
17 1/4"	0°, 4°	24" HP + 1 1/2" Spacer
	8°, 12°	24" + 1 1/2" Spacer
	0°, 4°	26"
17 1/2"	8°, 12°	26" + 3/4" Spacer
17 3/4"	0°, 4°	26" + 1 1/2" Spacer
18"	0°, 4°	26" + 3/4" Spacer
18 1/4"	8°, 12°	26" + 1 1/2" Spacer
18 3/4"	0°, 4°	26" + 1 1/2" Spacer

Front Seat Height	Caster Size	Fork & Stem Combination
17 1/4"	3"	5 1/4" + Std Stem
18"	3"	5 1/4" + 3/4" Stem
		6" + Std Stem
18 1/4"	5"	5 1/4" + Std Stem
18 3/4"	3"	5 1/4" + 1 1/2" Stem
		6" + 3/4" Stem
	6"	5 1/4" + Std Stem
19"	3"	7" + Std Stem
	5"	5 1/4" + 3/4" Stem
		6" + Std Stem
19 1/2"	3"	6" + 1 1/2" Stem
	6"	5 1/4" + 3/4" Stem
		6" + Std Stem
19 3/4"	3"	7" + 3/4" Stem
	5"	5 1/4" + 1 1/2" Stem
		6" + 3/4" Stem
20"	5"	7" + Std Stem
20 1/4"	6"	5 1/4" + 1 1/2" Stem
		6" + 3/4" Stem
20 1/2"	3"	7" + 1 1/2" Stem
	5"	6" + 1 1/2" Stem
	6"	7" + Std Stem
20 3/4"	5"	7" + 3/4" Stem
21"	6"	6" + 1 1/2" Stem

ACCESSORIES			
Backpack			
<input type="checkbox"/>	033A1	Black Utility Backpack	\$55.00
<input type="checkbox"/>	033A2	Black/Charcoal Kids	\$50.00
Seat Pouch			
<input type="checkbox"/>	033A8	Black	\$45.00
Extra Leg Strap			
<input type="checkbox"/>	033A55	K0038 Black	\$30.00
Miscellaneous			
<input type="checkbox"/>	033A82	Caddy	\$45.00
<input type="checkbox"/>	033A83	Tool Kit	\$25.00
Positioning Belt			
<input type="checkbox"/>	033A13	E0978 2" Aircraft Buckle Padded	\$120.00
<input type="checkbox"/>	033A12	E0978 2" Aircraft Buckle	\$110.00
<input type="checkbox"/>	033A142	E0978 1.5" Auto Buckle	\$70.00
Spoke Guard			
<input type="checkbox"/>	033A24	K0065 Clear	\$115.00
<input type="checkbox"/>	033A23	K0065 Black	\$115.00
Wheel Chair Tray Table			
<input type="checkbox"/>	033A30	E0950 Small Tuffak 13"-14"	\$250.00
<input type="checkbox"/>	033A31	E0950 Medium Tuffak 15"-17"	\$250.00
<input type="checkbox"/>	033A32	E0950 Large Tuffak 18"-20"	\$250.00
Side Guard			
<i>NA w/ height adj armrests or folding aluminum side guards</i>			
<i>Must be removed to fold backrest down</i>			
<input type="checkbox"/>	033SG3	K0108 Plastic Regular	\$125.00
<input type="checkbox"/>	033SG4	K0108 Plastic Kids	\$125.00
TOTAL CHAIR PRICE			\$

JAY SEATING OPTIONS*	
Select Cushion	Select Back
<input type="checkbox"/> Basic (E2601/E2602)	<input type="checkbox"/> Basic (E2611)
<input type="checkbox"/> Soft Combi P (E2605, E2606)	<input type="checkbox"/> J2 (E2613)
<input type="checkbox"/> Xtreme (E2603)	<input type="checkbox"/> J2 Deep Contour (E2620)
<input type="checkbox"/> Care (E2607)	<input type="checkbox"/> J2 Plus (E2613/E2614)
<input type="checkbox"/> Duo (E2603)	<input type="checkbox"/> Care (E2615)
<input type="checkbox"/> Easy (E2607)	<input type="checkbox"/> J3 (See OF)
<input type="checkbox"/> Jay 2 (K0734/K0735)	
<input type="checkbox"/> J2 Recline (K0734/K0735)	Indicate Backrest Width
<input type="checkbox"/> J2 Deep Contour (K0734.K0735)	_____
<input type="checkbox"/> J3 (K0734 or K0736)	Indicate Exact Cushion Dimensions
	(w" x d") _____

If a modification is necessary to achieve your cushion or backrest dimensions, do you want customer service to notify you prior to placing this order?

Please Choose

- Send Quote for any Modification
- Send Quote only when a modification has a charge
 - Fax me # _____
 - Email me _____

* Please see Jay Price List for latest pricing information. Also, note that Jay cushions and backs may ship from different distribution centers.