



Quickie M6™
January 2012
 EIF1

Mark For:

Submitting for: Quote Order

Date: _____
 Dealer Acct #: _____
 Dealer: _____
 Dealer Contact: _____
 Dealer Address: _____
 Dealer City: _____ ST: _____ ZIP: _____
 Dealer Phone: () _____ Fax: () _____
 Confirmation Email: _____
 Confirm Via: Fax Email

PO#: _____

ADDITIONAL SHIPPING INFORMATION

Ship To: _____
 Attention: _____
 Address: _____
 Address: _____
 Ship To City: _____ ST: _____ ZIP: _____
 Ship To Phone: () _____ Fax: () _____

The **HCPCS CODES** provided are based on code verification through PDAC or our interpretation of the code definitions.
 The use of **HCPCS CODES** does not ensure coverage or payment. For coverage information, verify the appropriate payer's coverage policy.

✓	Option #	HCPCS	Definition	Price	✓	Option #	HCPCS	Definition	Price
---	----------	-------	------------	-------	---	----------	-------	------------	-------

1. FRAME STYLE & SETUP

Base Model			
<input type="checkbox"/>	115WM21	K0009	Quickie M6 \$1,950.00

Frame Color			
<input type="checkbox"/>	115C2		Black NC
<input type="checkbox"/>	115C10		Black Opal NC
<input type="checkbox"/>	115C38		Blue Opal NC
<input type="checkbox"/>	115C36		Candy Blue NC
<input type="checkbox"/>	115C35		Candy Red NC
<input type="checkbox"/>	115C128		Copper NC
<input type="checkbox"/>	115C6		Forest Green NC
<input type="checkbox"/>	115C127		Fuchsia NC
<input type="checkbox"/>	115C126		Glow NC
<input type="checkbox"/>	115C4		Midnight Purple NC
<input type="checkbox"/>	115C5		Pearl Pink NC
<input type="checkbox"/>	115C103		Rootbeer NC
<input type="checkbox"/>	115C114		Sparkle Silver NC
<input type="checkbox"/>	115C113		Titanium Color Paint NC
<input type="checkbox"/>	115C11		Toxic Green NC
<input type="checkbox"/>	115C8		Yellow NC

Custom Frame Color			
<input type="checkbox"/>	115C112		Desert Camouflage \$295.00
<input type="checkbox"/>	115C43		Mossy Oak Camouflage \$295.00
<input type="checkbox"/>	115C42		Stars & Stripes \$295.00
<input type="checkbox"/>	115C108		Zebra Choose above base color. Std w/ white unless above base paint is chosen. NA with Black. \$295.00

Frame Width			
<input type="checkbox"/>	115W13	E2201	22" \$450.00
<input type="checkbox"/>	115W14	E2202	24" \$450.00
<input type="checkbox"/>	115W15	E2202	26" \$475.00
<input type="checkbox"/>	115W16	K0108	28" \$475.00
<input type="checkbox"/>	115W17	K0108	30" \$475.00

Frame Length			
<input type="checkbox"/>	115FL9		Regular (18"-19" depth) NC
<input type="checkbox"/>	115FL18		Long (20"-22" depth) NC

Seat Height			
<input type="checkbox"/>	115HT65		17" Includes low adjustable axle plate and shorter fork. Only available with 5" caster \$300.00
<input type="checkbox"/>	115HT11		18.5" Only available w/ 5" casters NC
<input type="checkbox"/>	115HT14		20" NC

Seat Depth			
<input type="checkbox"/>	115SD11		18" NC
<input type="checkbox"/>	115SD12		19" NC
<input type="checkbox"/>	115SD13	E2203	20" (Long frame only) \$475.00
<input type="checkbox"/>	115SD14	E2203	21" (Long frame only) \$475.00
<input type="checkbox"/>	115SD33	E2204	22" (Long frame only) \$795.00

Cushion			
<input type="checkbox"/>	115CU31	E2602	2" Cushion \$120.00

Backrest			
<input type="checkbox"/>	115BT34		8° w/ Push Handles Adjustable (16-3/4" to 20-1/8") NC

Back Height Preset			
<input type="checkbox"/>	115B91		16-3/4" NC
<input type="checkbox"/>	115B96		17-7/8" NC
<input type="checkbox"/>	115B39		19" NC
<input type="checkbox"/>	115B98		20-1/8" NC

Armrest			
<input type="checkbox"/>	115AA33	E0973	Dual Post Height Adj - Desk 10" \$245.00
<input type="checkbox"/>	115AA34	E0973	Dual Post Height Adj - Full 14" \$245.00

✓	Option #	HCPCS	Definition	Price	✓	Option #	HCPCS	Definition	Price
---	----------	-------	------------	-------	---	----------	-------	------------	-------

2. FRONT SEAT HEIGHT SET UP

Casters				
<input type="checkbox"/>	115CW6		5" Soft Roll	\$85.00
<input type="checkbox"/>	115CW36	E2214 E2215	8" x 2" Pneumatic	\$145.00
<input type="checkbox"/>	115CW37	E2214 K0108	8" x 2" Pneumatic w/Airless Insert	\$225.00

Front Hangers				
<input type="checkbox"/>	115H1		60° Swing Away	NC
<input type="checkbox"/>	115H146		70° Heavy-Duty Ft. Mount Lift Off <i>Only available with 5" caster and ang adj footplates. NA w/ extension tubes</i>	\$95.00
<input type="checkbox"/>	115H132		80° Heavy-Duty Ft. Mount Lift Off <i>Only available with 5" caster and ang adj footplates. NA w/ extension tubes</i>	\$95.00
<input type="checkbox"/>	115H41	E0990	ELR - Heavy-Duty Reinforced <i>Only available with aluminum footplates</i>	\$595.00

Footplate				
<input type="checkbox"/>	115F1	E0951	Composite w/ Heel Loops	NC
<input type="checkbox"/>	115F25	E0951 K0040	Adult Angle Adj w/ Heel Loops	\$185.00
<input type="checkbox"/>	115F5		Aluminum	NC

Extension Tubes				
<input type="checkbox"/>	115E1		Short (13"-15.5")	NC
<input type="checkbox"/>	115E2		Medium (15.5"-18")	NC
<input type="checkbox"/>	115E5		Long (18"-20.5") NA w/5" caster	NC

Lower Leg Length Range Matrix

Hanger	Ext Tube	Footplate		
		Composite	Alum Angle Adj	Aluminum
60°	Short	11.50" - 13.25"	11.75" - 14.00"	11.25" - 13.50"
	Medium	13.75" - 16.00"	14.00" - 16.25"	13.50" - 14.00"
	Long	16.00" - 16.50"	16.00" - 17.75"	NA
70°	No Ext Tube	NA	9.75" - 13.50"	NA
80°	No Ext Tube	NA	10.75" - 14.75"	NA
ELR	Short	NA	NA	13.00" - 17.00"
	Medium	NA	NA	14.25" - 17.00"

All values are ± 0.25" and based on flat seat heights and do not account for the minimum 2" ground clearance or footplate thickness. Full range not available with all seat heights. ELR values based on 70° hanger angle

Seat Height Matrix

Caster	Front/Rear
5"	17" (Different Axle & Fork)
5"	18.5"
8"	20"

3. REAR SEAT HEIGHT SET UP

Adjustable Axle Plate				
<input type="checkbox"/>	115AL2		Standard (for 18.5", 20" Seat Ht.)	NC
<input type="checkbox"/>	115AL28		Low (only for 17" Seat Ht.)	NC

Axle				
<input type="checkbox"/>	115AX1		Quick Release	NC

Rear Anti-Tip				
<input type="checkbox"/>	115AT1	E0971	Rear Anti-Tip	\$135.00
<input type="checkbox"/>	115AT3		Omit Rear Anti-Tip	NC

Wheel Rim/Size				
<input type="checkbox"/>	115RW2		24" Spoke Wheel	NC

Rear Wheel Tire				
<input type="checkbox"/>	115RT1	E2211 E2212	Pneumatic Tire	NC
<input type="checkbox"/>	115RT2	E2211 E2213	Pneumatic Airless Insert Tire	\$80.00

Handrim				
<input type="checkbox"/>	115HR1		Aluminum Anodized	NC
<input type="checkbox"/>	115HR4		Plastic Coated	\$95.00
<input type="checkbox"/>	115HR7		Omit Handrim	NC

Wheel Lock				
<input type="checkbox"/>	115WL1		High Mount Push Wheel Lock	NC
<input type="checkbox"/>	115WL12		Do Not Mount Wheel Lock	NC
<input type="checkbox"/>	115WL10		Omit Wheel Lock	NC

Extension Handle				
<input type="checkbox"/>	115EH1		6" Pair	\$75.00

4. ACCESSORIES

Positioning Belt				
<input type="checkbox"/>	115A12	E0978	Belt 2" Aircraft Buckle	\$110.00
<input type="checkbox"/>	115A13	E0978	Belt 2" Aircraft Buckle Padded	\$120.00
<input type="checkbox"/>	115A142	E0978	Belt 1.5" Auto Buckle	\$70.00

Miscellaneous				
<input type="checkbox"/>	115A83		Tool Kit	\$45.00