



## Quickie GP/GPV EIGP70, Effective March 2011

*The HCPCS codes provided are based on code verification by the PDAC or our interpretation of the code definitions.  
Proper use of HCPCS codes does not ensure coverage or payment. For coverage information, verify the appropriate payor's coverage policy.*

### ACCOUNT INFORMATION

Account No: \_\_\_\_\_ Drop Ship/Ship To: \_\_\_\_\_  
 Date: \_\_\_\_\_ Name: \_\_\_\_\_  
 P.O. #: \_\_\_\_\_ Address: \_\_\_\_\_  
 Buyer: \_\_\_\_\_ City: \_\_\_\_\_ State: \_\_\_\_\_  
 Marked For: \_\_\_\_\_ Zip Code: \_\_\_\_\_ Phone Number: \_\_\_\_\_

#### QUICKIE GP/GPV

- |  |          |
|--|----------|
| <input type="checkbox"/> 024WM23 K0005 Quickie GP-70 <sup>o</sup>  | \$995.00 |
| <input type="checkbox"/> 024WM23 K0005 Quickie GPV-70 <sup>o</sup> | \$995.00 |
| <i>NA 12"-13"</i>  |          |

#### Frame Angle

- |  |    |
|--|----|
| <input type="checkbox"/> 024FA6 60 <sup>o</sup>                          | NC |
| <input type="checkbox"/> 024FA7 70 <sup>o</sup>                          | NC |
| <input type="checkbox"/> 024FA1 75 <sup>o</sup> (NA with rigid footrest) | NC |
| <input type="checkbox"/> 024FA8 80 <sup>o</sup> (NA with rigid footrest) | NC |

#### Frame Color

- |   |    |
|---|----|
| <input type="checkbox"/> 024C2 Black                  | NC |
| <input type="checkbox"/> 024C10 Black Opal            | NC |
| <input type="checkbox"/> 024C38 Blue Opal             | NC |
| <input type="checkbox"/> 024C36 Candy Blue            | NC |
| <input type="checkbox"/> 024C35 Candy Red             | NC |
| <input type="checkbox"/> 024C128 Copper               | NC |
| <input type="checkbox"/> 024C127 Fuchsia              | NC |
| <input type="checkbox"/> 024C126 Glow                 | NC |
| <input type="checkbox"/> 024C6 Forest Green           | NC |
| <input type="checkbox"/> 024C4 Midnight Purple        | NC |
| <input type="checkbox"/> 024C5 Pearl Pink             | NC |
| <input type="checkbox"/> 024C103 Rootbeer             | NC |
| <input type="checkbox"/> 024C114 Sparkle Silver       | NC |
| <input type="checkbox"/> 024C113 Titanium Color Paint | NC |
| <input type="checkbox"/> 024C11 Toxic Green           | NC |
| <input type="checkbox"/> 024C8 Yellow                 | NC |

#### Custom Colors

- |  |          |
|--|----------|
| <input type="checkbox"/> 024C112 Desert Camouflage               | \$295.00 |
| <input type="checkbox"/> 024C42 Stars & Stripes                  | \$295.00 |
| <input type="checkbox"/> 024C43 Mossy Oak Camouflage             | \$295.00 |
| <input type="checkbox"/> 024C108 Zebra (choose base color above) | \$295.00 |

*Std w/ white unless paint color is chosen, NA w/ Black*

#### Xtender Power Assist

*Rear Seat Height of 16.75"*

- |  |            |
|--|------------|
| <input type="checkbox"/> 024XT1 E0986 Xtender Power Assist | \$6,595.00 |
|--|------------|

#### Xtender Pwr Asst Handrim

*Choose one of the following torque options*

- |  |    |
|--|----|
| <input type="checkbox"/> 024XH1 3x Torque Multiplier (Plastic Coated)    | NC |
| <input type="checkbox"/> 024XH2 1.5x Torque Multiplier (Stainless Steel) | NC |

#### Xtender Pwr Asst Tire

- |  |    |
|--|----|
| <input type="checkbox"/> 024XW1 E2211 Xtender Pneumatic Tire | NC |
| E2212  |    |

#### Xtender Pwr Asst Acc

- |   |            |
|---|------------|
| <input type="checkbox"/> 024XA1 Xtender Spare NiMH Battery Pack | \$1,395.00 |
|---|------------|

### FRAME DIMENSIONS

*Seat width 1/2" narrower*

- |                                    |    |
|------------------------------------|----|
| <input type="checkbox"/> 024W3 12" | NC |
| <input type="checkbox"/> 024W4 13" | NC |
| <input type="checkbox"/> 024W5 14" | NC |
| <input type="checkbox"/> 024W6 15" | NC |
| <input type="checkbox"/> 024W7 16" | NC |

#### Frame Width (Continued)

- |   |          |
|---|----------|
| <input type="checkbox"/> 024W8 17"        | NC       |
| <input type="checkbox"/> 024W9 18"        | NC       |
| <input type="checkbox"/> 024W41 19"       | \$125.00 |
| <input type="checkbox"/> 024W11 E2201 20" | \$450.00 |
| <input type="checkbox"/> 024W12 E2201 21" | \$450.00 |
| <input type="checkbox"/> 024W13 E2201 22" | \$475.00 |

#### Sling Depth

- |  |          |
|--|----------|
| <input type="checkbox"/> 024SD5 12"        | NC       |
| <input type="checkbox"/> 024SD6 13"        | NC       |
| <input type="checkbox"/> 024SD7 14"        | NC       |
| <input type="checkbox"/> 024SD8 15"        | NC       |
| <input type="checkbox"/> 024SD9 16"        | NC       |
| <input type="checkbox"/> 024SD10 17"       | NC       |
| <input type="checkbox"/> 024SD11 18"       | NC       |
| <input type="checkbox"/> 024SD72 19"       | \$75.00  |
| <input type="checkbox"/> 024SD13 E2203 20" | \$475.00 |
| <input type="checkbox"/> 024SD14 E2203 21" | \$475.00 |
| <input type="checkbox"/> 024SD33 E2204 22" | \$795.00 |

#### Frame Length

- |  |    |
|--|----|
| <input type="checkbox"/> 024FL30 18-1/4" (Avail. on 12" Sling Depth)           | NC |
| <input type="checkbox"/> 024FL31 19-1/4" (Avail. on 13" Sling Depth)           | NC |
| <input type="checkbox"/> 024FL32 20-1/4" (Avail. on 14", 15" Sling Depth)      | NC |
| <input type="checkbox"/> 024FL33 21-1/4" (Avail. on 14", 15", 16" Sling Depth) | NC |
| <input type="checkbox"/> 024FL34 22-1/4" (Avail. on 16", 17" Sling Depth)      | NC |
| <input type="checkbox"/> 024FL35 23-1/4" (Avail. on 17", 18" Sling Depth)      | NC |
| <input type="checkbox"/> 024FL36 24-1/4" (Avail. on 18" Sling Depth)           | NC |
| <input type="checkbox"/> 024FL37 25-1/4" (Avail. on 19" Sling Depth)           | NC |
| <input type="checkbox"/> 024FL38 26-1/4" (Avail. on 20" Sling Depth)           | NC |
| <input type="checkbox"/> 024FL39 27-1/4" (Avail. on 21" Sling Depth)           | NC |
| <input type="checkbox"/> 024FL40 28-1/4" (Avail. on 22" Sling Depth)           | NC |

#### Frame Height (See Chart)

- |  |          |
|--|----------|
| <input type="checkbox"/> 024FH2 Short (front 9.25"-rear 8.25") | \$110.00 |
| <i>NA with Xtender Power Assist</i>                            |          |
| <input type="checkbox"/> 024FH5 Standard (front 10" - rear 9") | NC       |
| <input type="checkbox"/> 024FH7 Tall (front 11" - rear 10")    | \$110.00 |

### SEATING

#### Back Type

- |  |                                     |    |
|--|-------------------------------------|----|
| <input type="checkbox"/> 024BT3                      | <b>Non-Folding Backrest (Fixed)</b> |    |
| <i>NA with push handles or adjustable upholstery</i> |                                     |    |
| <input type="checkbox"/> 024BB1 8" Back Height       |                                     | NC |
| <input type="checkbox"/> 024BB2 9" Back Height       |                                     | NC |
| <input type="checkbox"/> 024BB3 10" Back Height      |                                     | NC |
| <input type="checkbox"/> 024BB4 11" Back Height      |                                     | NC |
| <input type="checkbox"/> 024BB5 12" Back Height      |                                     | NC |
| <input type="checkbox"/> 024BB6 13" Back Height      |                                     | NC |
| <input type="checkbox"/> 024BB7 14" Back Height      |                                     | NC |
| <input type="checkbox"/> 024BB8 15" Back Height      |                                     | NC |
| <input type="checkbox"/> 024BB9 16" Back Height      |                                     | NC |

**Back Type (Continued)**

<input type="checkbox"/> 024BT1	<b>Folding Backrest (Adjustable)</b>	
<input type="checkbox"/> 024BA1	11" - 14" Back Height	\$65.00
<input type="checkbox"/> 024BA2	12" - 16" Back Height	\$65.00
<input type="checkbox"/> 024BA3	14" - 18" Back Height	\$65.00
<input type="checkbox"/> 024BA4	16" - 20" Back Height	\$65.00

**Back Upholstery**

<input type="checkbox"/> 024BU1	Standard	NC
<input type="checkbox"/> 024BU8	Adjustable	\$175.00
<input type="checkbox"/> 024BU16	Reinforced (NA with adjustable upholstery)	\$30.00
<input type="checkbox"/> 024BU14	Omit Back Upholstery	-\$15.00

**Backrest Options**

<input type="checkbox"/> 024BT33	Folding w/ Integral Push Handles	NC
<input type="checkbox"/> 024A85	Bolt-On Push Handle (Folding only) <i>NA with Xtender Power Assist</i>	\$115.00
<input type="checkbox"/> 024SR2	Regular Stroller Handle <i>NA with amputee axle plate or Xtender Power Assist</i>	\$125.00
<input type="checkbox"/> 024SR3	Tall Stroller Handle <i>NA with amputee axle plate or Xtender Power Assist</i>	\$125.00

**Cushion**

<input type="checkbox"/> 024CU9	Reinforced Seat Sling	\$30.00
<input type="checkbox"/> 024CU1	K0669 2" Cushion	\$95.00
<input type="checkbox"/> 024CU2	Omit Cushion	NC
<input type="checkbox"/> 024CU4	Omit Cushion and Sling	-\$15.00

**ARMREST**

<input type="checkbox"/> 024AA9	E0973 <b>Single Post Ht. Adjust. Low</b>	\$260.00
<input type="checkbox"/>	<i>Must be removed to fold back down</i>	
<input type="checkbox"/> 024AM1	E0973 Desk Length (10")	
<input type="checkbox"/> 024AM2	E0973 Full Length (14")	
<input type="checkbox"/> 024AA8	E0973 <b>Single Post Ht. Adjust. Std</b>	\$260.00
<input type="checkbox"/>	<i>Must be removed to fold back down</i>	
<input type="checkbox"/> 024AM1	E0973 Desk Length (10")	
<input type="checkbox"/> 024AM2	E0973 Full Length (14")	
<input type="checkbox"/>	E0973 <b>Padded Swing-Away Armrest</b>	
	<i>All Padded Swing-Away NA with amputee axle plate or Xtender Power Assist</i>	
<input type="checkbox"/> 024AA19	E0973 Padded S-A Medium	\$175.00
<input type="checkbox"/> 024AA16	E0973 Padded S-A Adult	\$175.00

**FOOTREST**

**1. GP Footrest Options**

<input type="checkbox"/> 024F56	Rigid w/ Leg Strap (NA w/ Ext. Tubes)	NC
	<i>NA with 75° or 80° frame angle</i>	
<input type="checkbox"/> 024F50	K0038 Platform w/ Leg Strap (NA w/ Ext. Tubes)	\$95.00
<input type="checkbox"/> 024F1	Composite Footplate w/ Heel Loops	NC
<input type="checkbox"/> 024F3	K0040 Adult Angle Adj. w/ Heel Loops	\$185.00
	E0951	
<input type="checkbox"/> 024F48	K0038 Platform Flip Up w/Leg Strap	\$190.00

**2. GPV Footrest Options**

<input type="checkbox"/> 024F56	Rigid w/Leg Strap (NA w/ Ext. Tubes)	NC
	<i>NA with 75° or 80° frame angle</i>	
<input type="checkbox"/> 024F50	K0040 Platform w/ Leg Strap (NA w/ Ext. Tubes)	\$95.00
<input type="checkbox"/> 024F48	K0038 Platform Flip Up w/Leg Strap	\$190.00
	K0040	

**Extension Tubes**

<input type="checkbox"/> 024E1	Short (14" to 16-1/2")	NC
<input type="checkbox"/> 024E2	Medium (16-1/2" to 19")	NC
<input type="checkbox"/> 024E5	Long (19" to 21-1/2")	NC
<input type="checkbox"/> 024E16	K0037 High Mount	\$55.00
	<i>Only available with adult angle adjustable footplates</i>	

**Footrest Options**

<input type="checkbox"/> 024IG2	Impact Guard Neoprene	\$85.00
<input type="checkbox"/> 024HL4	Omit Heel Loop	NC
<input type="checkbox"/> 024HL6	Omit Single Leg Strap	NC
<input type="checkbox"/> 024F7	Omit Footplate & Extension Tube	-\$15.00

**CASTERS**

**Casters**

<input type="checkbox"/> 024CW1	3" Micro	NC
<input type="checkbox"/> 024CW3	4" Poly	\$65.00
<input type="checkbox"/> 024CW31	4" Micro Lighted	\$70.00
<input type="checkbox"/> 024CW6	5" Soft Roll Semi-Pneumatic	\$85.00
<input type="checkbox"/> 024CW5	5" Low Profile Poly	NC
<input type="checkbox"/> 024CW7	5" x 1.5" Semi-Pneumatic	\$85.00
	<i>Only available for use in bottom fork holes</i>	
<input type="checkbox"/> 024CW33	E2214 6" Pneumatic	\$75.00
	E2215	
<input type="checkbox"/> 024CW34	6" Polyurethane	NC
<input type="checkbox"/> 024CW35	6" x 1.5" Semi-Pneumatic	\$85.00
	<i>Only available for use in bottom fork holes</i>	
<input type="checkbox"/> 024CW21	E2214 8" Pneumatic	\$75.00
	E2215 <i>Only available on 60° frame option</i>	
<input type="checkbox"/> 024CW32	8" Polyurethane	NC
	<i>Only available on 60° frame option</i>	
<input type="checkbox"/> 024CW36	E2214 8" x 2" Pneumatic	\$145.00
	E2215 <i>Only available on 60° frame option</i>	
<input type="checkbox"/> 024CW37	E2214 8" x 2" Pneumatic w/ airless inserts	\$225.00
	K0108 <i>Only available on 60° frame option</i>	

**Caster Options**

<input type="checkbox"/> 024CS5	3/4" Fork Stem	NC
<input type="checkbox"/> 024CS9	1-1/2" Fork Stem	NC
<input type="checkbox"/> 024CS3	K0073 Caster Pin Locks	\$195.00
<input type="checkbox"/> 024CW12	Omit Caster Wheel & Tire	NC

**AXLE AND WHEEL**

**Axle Plate**

*NA with Xtender Power Assist*

<input type="checkbox"/> 024AL2	Standard Adjustable	NC
<input type="checkbox"/> 024AL13	Double Slotted Adjustable	\$60.00
<input type="checkbox"/> 024AL6	E0959 Amputee Adjustable	\$175.00
	<i>NA with padded swing-away armrests or stroller handles</i>	

**Axle**

<input type="checkbox"/> 024AX1	Quick Release	NC
<input type="checkbox"/> 024AX11	Quad Quick Release	\$90.00

**Anti Tips**

<input type="checkbox"/> 024AT1	E0971 Rear Anti-Tip	\$135.00
---------------------------------	---------------------	----------

**Wheels**

<input type="checkbox"/> 024RW8	22" Mag (NA with Xtender Power Assist)	NC
<input type="checkbox"/> 024RW9	22" Spoke (NA with Xtender Power Assist)	\$200.00
<input type="checkbox"/> 024RW1	24" Mag	NC
<input type="checkbox"/> 024RW2	24" Spoke	NC
<input type="checkbox"/> 024RW11	26" Spoke (NA with Xtender Power Assist)	\$65.00
<input type="checkbox"/> 024RR25	Omit Rear Wheel and Axle	-\$55.00
	(Intended size for set up _____)	

**Quickie Performance**

*NA with Xtender Power Assist, spoke guards or titanium axles*

*Contact Built 4 Me for Spinergy color options*

<input type="checkbox"/> 024RW42	24" Spinergy (540)(Clinchers)	\$695.00
	<i>Available w/ Black Spokes only</i>	

**Tire**

<input type="checkbox"/> 024RT1	E2211 Pneumatic Tire	NC
	E2212 <i>NA with heat-treated wheels</i>	
<input type="checkbox"/> 024RT2	E2213 Pneumatic Airless Insert Tire	\$80.00
	E2211 <i>NA with heat-treated wheels or 26" wheels</i>	
<input type="checkbox"/> 024RT3	E2218 Full Polyurethane Tire	NC
	<i>NA with heat-treated wheels or 26" wheels</i>	
<input type="checkbox"/> 024RT6	E2211 High-Pressure Clincher (24"/26" Only)	\$95.00
	E2212 <i>NA with heat-treated wheels or mag wheels</i>	

Handrim		
<input type="checkbox"/>	024HR1	Aluminum Anodized NC
<input type="checkbox"/>	024HR4	Plastic Coated \$95.00
<input type="checkbox"/>	024HR7	Omit Handrim NC

Projections		
<i>NA on Plastic Coated Handrims</i>		
<input type="checkbox"/>	024HR8	8 Oblique Handrim (24"/26") \$240.00

Wheel Locks		
<input type="checkbox"/>	024WL1	High Mount Push NC
<input type="checkbox"/>	024WL2	High Mount Pull NC
<input type="checkbox"/>	024WL5	High Mount Scissor NC
<i>NA with grade aids or extension handles</i>		
<input type="checkbox"/>	024WL12	Do Not Mount Wheel Lock NC
<input type="checkbox"/>	024WL10	Omit Wheel Lock NC

Extension Handle		
<input type="checkbox"/>	024EH1	E0961 6" Pair \$75.00

Grade Aid		
<i>NA with poly, turbo, high pressure tires or Xtender Power Assist</i>		
<input type="checkbox"/>	024GA1	E0974 Grade Aid \$150.00

ACCESSORIES		
Back Pack		
<input type="checkbox"/>	024A1	Black Utility Backpack \$55.00
<input type="checkbox"/>	024A2	Black/Charcoal Kids \$50.00

Seat Pouch		
<input type="checkbox"/>	024A8	Black \$45.00

Crutch Holder		
<input type="checkbox"/>	024A53	E2207 Left Side \$195.00
<input type="checkbox"/>	024A52	E2207 Right Side \$195.00

Extra Leg Strap		
<input type="checkbox"/>	024A55	K0038 Extra Leg Strap Black \$30.00

Miscellaneous		
<input type="checkbox"/>	024A83	Tool Kit \$45.00

Positioning Belt		
<input type="checkbox"/>	024A142	E0978 1.5" Auto Buckle \$70.00
<input type="checkbox"/>	024A12	E0978 2" Aircraft Buckle \$110.00
<input type="checkbox"/>	024A13	E0978 2" Aircraft Buckle Padded \$120.00

Spoke Guard		
<i>NA with Spineray Wheels or Xtender Power Assist</i>		
<input type="checkbox"/>	024A23	K0065 Black \$115.00
<input type="checkbox"/>	024A24	K0065 Clear \$115.00

Tray Table		
<input type="checkbox"/>	024A31	E0950 Tray Table \$250.00

Side Guard		
<input type="checkbox"/>	024SG3	K0108 Plastic Regular \$125.00
<input type="checkbox"/>	024SG4	K1008 Plastic Kids \$125.00

JAY SEATING OPTIONS	
Select Cushion	Select Back
<input type="checkbox"/> Basic (E2601/E2602)	<input type="checkbox"/> Basic (E2611)
<input type="checkbox"/> Soft Combi P (E2605, E2606)	<input type="checkbox"/> J2 (E2613)
<input type="checkbox"/> Xtreme (E2603)	<input type="checkbox"/> J2 Deep Contour (E2620)
<input type="checkbox"/> Care (E2607)	<input type="checkbox"/> Xtreme (K0669)
<input type="checkbox"/> Duo (E2603)	<input type="checkbox"/> J2 Plus (E2613/E2614)
<input type="checkbox"/> Easy (E2607)	<input type="checkbox"/> Care (E2615)
<input type="checkbox"/> Jay 2 (K0734)	<input type="checkbox"/> J3 (See OF)
<input type="checkbox"/> J2 Recline (K0734)	
<input type="checkbox"/> J2 Deep Contour (K0734)	
<input type="checkbox"/> J3 Cushion (K0734 or K0736)	

Indicate Backrest Width \_\_\_\_\_

Indicate Exact Cushion Dimensions

(w x d) \_\_\_\_\_

If a modification is necessary to achieve your cushion or backrest dimensions, do you want customer service to notify you prior to placing this order?

**Please Choose**

- Send Quote for any modification
- Send Quote only when a modification has a charge


Fax me # \_\_\_\_\_

Email me \_\_\_\_\_

\* Please see Jay Price List for latest pricing information.

**TOTAL CHAIR PRICE** \$ \_\_\_\_\_

**Frame Dimensions**



### GP Series Front STF Chart

*\*All measurements given are an estimate and have been rounded to the nearest .5"*

SHORT Frame ( 9 1/4" Front / 8 1/4" Rear)										
Fork	Spacer	Fork Hole	3"	4" Poly	5" Poly	6" Poly	6" Pneu	8" Poly	8" Pneu	8x2
				4" Lighted	5" x 1.5					
4 1/4	STD	Bottom	15.50	16.00	16.50	17.00	17.00	x	x	x
4 1/4	3/4	Bottom	16.50	17.00	17.00	17.50	18.00	x	x	x
4 1/4	1 1/2	Bottom	17.00	17.50	18.00	18.50	18.50	x	x	x
5 1/4	STD	Top	x	x	16.50	17.00	17.00	x	x	x
5 1/4	STD	Bottom	16.50	17.00	17.5	18.00	18.00	19.00	20.5	x
5 1/4	3/4	Top	x	x	17.50	18.00	18.00	x	x	x
5 1/4	3/4	Bottom	17.50	18.00	18.50	18.50	19.00	19.50	20.00	x
5 1/4	1 1/2	Top	x	x	18.00	18.50	18.50	x	x	x
5 1/4	1 1/2	Bottom	18.00	18.50	19.00	19.50	19.50	20.50	20.50	x
6	STD	Top	x	x	17.50	18.00	18.00	18.50	19.00	x
6	STD	Bottom	17.50	18.00	18.50	19.00	19.00	19.50	20.00	x
6	3/4	Top	x	x	18.00	18.50	18.50	19.50	19.50	x
6	3/4	Bottom	18.00	18.50	19.00	19.50	19.50	20.50	20.50	x
6	1 1/2	Top	x	x	19.00	19.50	19.50	20.00	20.50	x
6	1 1/2	Bottom	19.00	19.50	20.00	20.50	20.50	21.00	21.50	x
7	STD	Top	x	x	18.50	19.00	19.00	19.50	20.00	x
7	STD	Bottom	18.50	19.00	19.50	20.00	20.00	20.50	21.00	x
7	3/4	Top	x	x	19.00	19.50	19.50	20.50	20.50	x
7	3/4	Bottom	19.00	19.50	20.00	20.50	20.50	21.50	21.50	x
7	1 1/2	Top	x	x	20.00	20.50	20.50	21.00	21.50	x
7	1 1/2	Bottom	20.00	20.50	21.00	21.00	21.50	22.00	22.50	x
8x2	STD	Top	x	x	x	x	x	x	x	19.00
8x2	STD	Bottom	x	x	x	x	x	x	x	20.00
8x2	3/4	Top	x	x	x	x	x	x	x	19.50
8x2	3/4	Bottom	x	x	x	x	x	x	x	20.50
8x2	1 1/2	Top	x	x	x	x	x	x	x	20.50
8x2	1 1/2	Bottom	x	x	x	x	x	x	x	21.00

<b>STANDARD Frame (10" Front / 9" Rear) ADD 3/4" to each measurement above</b>
<b>TALL Frame (11" Front / 10" Rear) ADD 1 3/4" to each measurement above</b>

\*\*\* Custom Frame Height through Configurator - Front 9 - 13" - Calculate Accordingly (Minimum Front STF = 15 1/2" / Max = 26")

### GP Series Rear STF Chart

*\*All measurements given are an estimate and have been rounded to the nearest .5"*

SHORT Frame ( 9 1/4" Front / 8 1/4" Rear)			Rear Axle Plate Hole Position					
Wheel	Rim	Tire	Standard			Inverted		
			Top	Middle	Bottom	Top	Middle	Bottom
22"	Mag / Spoke	Full Poly	14.5	15	15.5	16	16.5	17
22"	Mag / Spoke	Pneum / Insert	14.5	15	16	16	17	17.5
24"	Mag / Spoke	Full Poly	15	15.5	16.5	16.5	17	17.5
24"	Spoke / Spinerger	Clincher	15	16	16.5	16.5	17	18
24"	Mag / Spoke / Spinerger	Pneum / Insert	15.5	16	17	17	17.5	18
26"	Spoke / Spinerger	Clincher	16	17	17.5	17.5	18	19
26"	Spoke / Spinerger	Pneumatic	16.5	17	17.5	18	18.5	19

TALL Frame ( 11" Front / 10" Rear)			Rear Axle Plate Hole Position					
Wheel	Rim	Tire	Standard			Inverted		
			Top	Middle	Bottom	Top	Middle	Bottom
22"	Mag / Spoke	Full Poly	16	16.5	17.5	17.5	18	19
22"	Mag / Spoke	Pneum / Insert	16.5	17	17.5	18	18.5	19
24"	Mag / Spoke	Full Poly	17	17.5	18	18.5	18.5	19.5
24"	Spoke / Spinerger	Clincher	17	17.5	18	18.5	19	19.5
24"	Mag / Spoke / Spinerger	Pneum / Insert	17	17.5	18.5	19	19	20
26"	Spoke / Spinerger	Clincher	18	18.5	19	19.5	20	20.5
26"	Spoke / Spinerger	Pneumatic	18	19	19.5	19.5	20	21

STANDARD Frame ( 10" Front / 9" Rear)			Rear Axle Plate Hole Position					
Wheel	Rim	Tire	Standard			Inverted		
			Top	Middle	Bottom	Top	Middle	Bottom
22"	Mag / Spoke	Full Poly	15	15.5	16.5	16.5	17.5	18
22"	Mag / Spoke	Pneum / Insert	15.5	16	17	17	17.5	18
24"	Mag / Spoke	Full Poly	16	16.5	17	17	18	18.5
24"	Spoke / Spinerger	Clincher	16	16.5	17	17.5	18	18.5
24"	Mag / Spoke / Spinerger	Pneum / Insert	16	16.5	17.5	17.5	18	18.5
26"	Spoke / Spinerger	Clincher	17	17.5	18	18.5	19	19.5
26"	Spoke / Spinerger	Pneumatic	17	18	18.5	18.5	19.5	20

Maximum Dump 2 3/4"  
 \*\*\* Custom Frame Height Through Configurator - REAR 6 1/2 - 13" (Max frame squeeze = 3 1/2") Calculate Accordingly