



# Quickie GP Swing-Away

EIGPSA, Effective March 2011

The HCPCS codes provided are based on code verification by the PDAC or our interpretation of the code definitions. Proper use of HCPCS codes does not ensure coverage or payment. For coverage information, verify the appropriate payor's coverage policy.

## ACCOUNT INFORMATION

Account No: _____	Drop Ship/Ship To: _____
Date: _____	Name: _____
P.O. #: _____	Address: _____
Buyer: _____	City: _____ State: _____
Marked For: _____	Zip Code: _____ Phone Number: _____

### QUICKIE GP SWING-AWAY

025WM26 K0005 Quickie GP Swing-Away \$1,327.00

### FRAMES

#### Frame Color

<input type="checkbox"/> 025C2	Black	NC
<input type="checkbox"/> 025C10	Black Opal	NC
<input type="checkbox"/> 025C38	Blue Opal	NC
<input type="checkbox"/> 025C36	Candy Blue	NC
<input type="checkbox"/> 025C35	Candy Red	NC
<input type="checkbox"/> 025C128	Copper	NC
<input type="checkbox"/> 025C6	Forest Green	NC
<input type="checkbox"/> 025C127	Fuchsia	NC
<input type="checkbox"/> 025C126	Glow	NC
<input type="checkbox"/> 025C4	Midnight Purple	NC
<input type="checkbox"/> 025C5	Pearl Pink	NC
<input type="checkbox"/> 025C103	Rootbeer	NC
<input type="checkbox"/> 025C114	Sparkle Silver	NC
<input type="checkbox"/> 025C113	Titanium Color Paint	NC
<input type="checkbox"/> 025C11	Toxic Green	NC
<input type="checkbox"/> 025C8	Yellow	NC

#### Custom Color

<input type="checkbox"/> 025C112	Desert Camouflage	\$295.00
<input type="checkbox"/> 025C42	Stars and Stripes	\$295.00
<input type="checkbox"/> 025C43	Mossy Oak Camouflage	\$295.00
<input type="checkbox"/> 025C108	Zebra (choose base color above)	\$295.00

*Std w/ white unless paint color is chosen, NA w/ Black*

#### Frame Width

*Seat width 1/2" narrower*

<input type="checkbox"/> 025W3	12" NA with adjustable upholstery	NC
<input type="checkbox"/> 025W4	13" NA with adjustable upholstery	NC
<input type="checkbox"/> 025W5	14"	NC
<input type="checkbox"/> 025W6	15"	NC
<input type="checkbox"/> 025W7	16"	NC
<input type="checkbox"/> 025W8	17"	NC
<input type="checkbox"/> 025W9	18"	NC
<input type="checkbox"/> 025W41	19"	\$125.00
<input type="checkbox"/> 025W11	E2201 20"	\$450.00

#### Sling Depth

<input type="checkbox"/> 025SD5	12"	NC
<input type="checkbox"/> 025SD6	13"	NC
<input type="checkbox"/> 025SD7	14"	NC
<input type="checkbox"/> 025SD8	15"	NC
<input type="checkbox"/> 025SD9	16"	NC
<input type="checkbox"/> 025SD10	17"	NC
<input type="checkbox"/> 025SD11	18"	NC
<input type="checkbox"/> 025SD74	19"	\$125.00
<input type="checkbox"/> 025SD13	E2203 20"	\$475.00
<input type="checkbox"/> 025SD14	E2203 21"	\$475.00
<input type="checkbox"/> 025SD33	E2204 22"	\$795.00

### SEATING

#### Cushion

<input type="checkbox"/> 025CU1	K0669 2"	\$95.00
<input type="checkbox"/> 025CU2	Omit Cushion	NC
<input type="checkbox"/> 025CU4	Omit Sling	-\$15.00

### BACKREST

#### Back Type

<input type="checkbox"/> 025BT1	Folding Backrest (Adjustable)	
<input type="checkbox"/> 025BA1	11"-14" Back Height	\$65.00
<input type="checkbox"/> 025BA2	12"-16" Back Height	\$65.00
<input type="checkbox"/> 025BA3	14"-18" Back Height	\$65.00
<input type="checkbox"/> 025BA4	16"-20" Back Height	\$65.00

#### Non-Folding Backrest (fixed)

*NA with adjustable upholstery or push handles*

<input type="checkbox"/> 025BB1	8" Back Height	NC
<input type="checkbox"/> 025BB2	9" Back Height	NC
<input type="checkbox"/> 025BB3	10" Back Height	NC
<input type="checkbox"/> 025BB4	11" Back Height	NC
<input type="checkbox"/> 025BB5	12" Back Height	NC
<input type="checkbox"/> 025BB6	13" Back Height	NC
<input type="checkbox"/> 025BB7	14" Back Height	NC
<input type="checkbox"/> 025BB8	15" Back Height	NC
<input type="checkbox"/> 025BB9	16" Back Height	NC

#### Back Upholstery

<input type="checkbox"/> 025BU1	Standard	NC
<input type="checkbox"/> 025BU8	Adjustable	\$175.00
	<i>NA with fabric side guards</i>	
<input type="checkbox"/> 025BU16	Reinforced	\$30.00
	<i>NA with adjustable upholstery</i>	
<input type="checkbox"/> 025BU14	Omit Upholstery	-\$15.00

#### Backrest Options

<input type="checkbox"/> 025BT33	Folding w/ Integral Push Handles	NC
<input type="checkbox"/> 025A85	Bolt-On Push Handles	\$115.00
	<i>NA with Xtender Power Assist</i>	
<input type="checkbox"/> 025SR2	Regular Stroller Handle	\$125.00
	<i>NA with Xtender Power Assist</i>	
<input type="checkbox"/> 025SR3	Tail Stroller Handle	\$125.00
	<i>NA with amputee axle plate or Xtender Power Assist</i>	

**CASTERS**

<input type="checkbox"/>	025CW1	3" Micro	NC
<input type="checkbox"/>	025CW3	4" Low Profile Polyurethane	\$65.00
<input type="checkbox"/>	025CW31	4" Micro Lighted	\$70.00
<input type="checkbox"/>	025CW5	5" Low Profile Polyurethane	NC
<input type="checkbox"/>	025CW7	E2219 5" x 1.5" Semi-Pneumatic	\$85.00
		<i>Only available for use in bottom fork hole</i>	
<input type="checkbox"/>	025CW6	5" Soft Roll Semi-Pneumatic	\$85.00
<input type="checkbox"/>	025CW33	E2214 6" Pneumatic	\$75.00
<input type="checkbox"/>	025CW34	6" Polyurethane	NC
<input type="checkbox"/>	025CW35	E2219 6" x 1.5" Semi-Pneumatic	\$85.00
		<i>Only available for use in bottom fork hole</i>	
<input type="checkbox"/>	025CW21	E2214 8" Pneumatic	\$75.00
<input type="checkbox"/>	025CW32	8" Polyurethane	NC
<input type="checkbox"/>	025CW36	E2214 8" x 2" Pneumatic	\$145.00
<input type="checkbox"/>	025CW37	E2214 8" x 2" Pneumatic w/Airless Inserts	\$225.00
		K0108	

**Caster Options**

*NA with Xtender Power Assist*

<input type="checkbox"/>	025CS5	3/4" Longer Stem Bolt	NC
<input type="checkbox"/>	025CS9	1-1/2" Longer Stem Bolt	NC
<input type="checkbox"/>	025OC1	K0073 Caster Pin Locks	\$195.00
		<i>NA with Xtender Power Assist</i>	

**ARMREST**

<input type="checkbox"/>	025AA9	E0973 <b>Single Post Ht. Adjust. Low</b>	\$260.00
<input type="checkbox"/>	025AM1	E0973 Desk Length (10")	
<input type="checkbox"/>	025AM2	E0973 Full Length (14")	
<input type="checkbox"/>	025AA8	E0973 <b>Single Post Ht. Adjust. Std</b>	\$260.00
<input type="checkbox"/>	025AM1	E0973 Desk Length (10")	
<input type="checkbox"/>	025AM2	E0973 Full Length (14")	
<input type="checkbox"/>		E0973 <b>Padded Swing-Away Armrest</b>	
		<i>All Padded Swing Away -NA with amputee axle plate or Xtender Power Assist</i>	
<input type="checkbox"/>	025AA19	E0973 Padded S-A Medium	\$175.00
<input type="checkbox"/>	025AA16	E0973 Padded S-A Adult	\$175.00

**FOOTREST**

**Hanger**

<input type="checkbox"/>	025H1	60° Swing Away	NC
<input type="checkbox"/>	025H2	70° Swing Away	NC
<input type="checkbox"/>	025H6	90° Swing Away	NC
		<i>Only available with kids angle adjustable, adult angle adj or 90° footboard footplates. NA with extension tubes</i>	
<input type="checkbox"/>	025H24	K0053 Articulating	\$350.00
		<i>Only available with composite and angle adjustable footplates</i>	
<input type="checkbox"/>	025H14	E0990 Elevating	\$295.00
		<i>Only available with aluminum or adult adjustable footplates</i>	
<input type="checkbox"/>	025H32	Omit	-\$15.00

**Footplate**

<input type="checkbox"/>	025F1	Composite w/ Heel Loops 14"-20"	NC
		<i>NA with 90° hangers and elevating legrest</i>	
<input type="checkbox"/>	025F3	K0040 Adult Angle Adj. w/ Heel Loops	\$185.00
		E0951	
<input type="checkbox"/>	025F4	K0040 Kids Angle Adj. w/ Leg Strap	\$175.00
		K0038 <i>Only available on 12" - 16" width frames.</i>	
<input type="checkbox"/>	025F5	Aluminum	NC
		<i>Only available with elevating hangers</i>	
<input type="checkbox"/>	025F48	K0038 Platform Flip Up w/Leg Strap	\$190.00
		K0040 <i>Only available on 60° and 70° hangers</i>	
<input type="checkbox"/>	025F7	Omit Footplate & Extension Tube	-\$15.00

**Extension Tube**

<input type="checkbox"/>	025E1	Short (See chart for lengths)	NC
		<i>NA with elevating or articulating legrests</i>	
<input type="checkbox"/>	025E2	Medium (See chart for lengths)	NC
<input type="checkbox"/>	025E5	Long (See chart for lengths)	NC
<input type="checkbox"/>	025E12	Short ELR	NC
<input type="checkbox"/>	025E15	Long ELR	NC
<input type="checkbox"/>	025E16	K0037 High Mount	\$55.00
		<i>NA with elevating and 90° elevating legrests</i>	
		<i>Only available with adult or kids angle adj footplates</i>	

**Footrest Accessories**

<input type="checkbox"/>	025HL4	Omit Heel Loop	NC
--------------------------	--------	----------------	----

**Impact Guard**

<input type="checkbox"/>	025IG2	Neoprene	\$85.00
		<i>Only available on 60° or 70° hangers</i>	

**AXLE AND WHEEL**

**Axle Plate**

*NA with 8" x 2" casters or Xtender Power Assist*

<input type="checkbox"/>	025AL2	Standard Adjustable	NC
<input type="checkbox"/>	025AL6	E0959 Amputee Adjustable	\$175.00
		<i>NA with padded swing-away armrests or stroller handles</i>	
<input type="checkbox"/>	025AL13	Double Slotted Adjustable	\$60.00

**Axle**

*NA with 8" x 2" casters or Xtender Power Assist*

<input type="checkbox"/>	025AX4	Quick Release Steel	NC
<input type="checkbox"/>	025AX11	Quad Quick Release	\$90.00

**Rear Anti-Tip**

<input type="checkbox"/>	025AT1	E0971 Rear Anti-Tip	\$135.00
--------------------------	--------	---------------------	----------

**Rear Wheel**

<input type="checkbox"/>	025RW8	22" Mag	NC
		<i>NA with 8" x 2" casters or Xtender Power Assist</i>	
<input type="checkbox"/>	025RW9	22" Spoke	\$200.00
		<i>NA with 8" x 2" casters or Xtender Power Assist</i>	
<input type="checkbox"/>	025RW1	24" Mag	NC
<input type="checkbox"/>	025RW2	24" Spoke	NC
<input type="checkbox"/>	025RW11	26" Spoke	\$65.00
		<i>NA with 8" x 2" casters or Xtender Power Assist</i>	

**Seat Sling to Top of Footplate Range**

60° Swing- Away	Short Ext.	12 - 15 <sup>3/4</sup>
60° Swing- Away	Medium Ext.	14 <sup>1/2</sup> - 18 <sup>1/4</sup>
60° Swing- Away	Long Ext.	17 - 20 <sup>3/4</sup>
60° Platform Flip-Up	Short Ext.	12 <sup>3/8</sup> - 15 <sup>5/8</sup>
60° Platform Flip-Up	Medium Ext.	14 <sup>1/2</sup> - 17 <sup>3/4</sup>
60° Platform Flip-Up	Long Ext.	17 - 20 <sup>1/4</sup>
70° Swing- Away	Short Ext.	11 - 14 <sup>3/4</sup>
70° Swing- Away	Medium Ext.	13 <sup>1/2</sup> - 17 <sup>1/4</sup>
70° Swing- Away	Long Ext.	16 - 19 <sup>3/4</sup>
70° Platform Flip-Up	Short Ext.	11 <sup>5/8</sup> - 15 <sup>3/8</sup>
70° Platform Flip-Up	Medium Ext.	13 <sup>3/4</sup> - 17 <sup>1/2</sup>
70° Platform Flip-Up	Long Ext.	16 <sup>1/4</sup> - 20

\*All adjustment ranges may vary due to seat height restrictions. Recommended 2<sup>1/2</sup>" footrest clearance to ground surface.

**Rear Wheel Tire**

- 025RT1 E2211 Pneumatic NC  
E2212
- 025RT2 E2213 Pneumatic Airless Insert \$80.00  
E2211 *NA on 26" wheels*
- 025RT3 E2218 Full Polyurethane NC  
*NA on 26" wheels*
- 025RT6 E2211 High-Pressure Clincher \$95.00  
E2212 *Only available on 24" and 26" spoke rims*  
*Not available for heat treated rims*
- 025RT5 E2211 Turbo Tread NC  
E2212 *Only available with 24" or 26" heat treated wheel rims*

**Handrim**

*NA with Xtender Power Assist*

- 025HR1 Aluminum Anodized NC
- 025HR4 Plastic Coated \$95.00
- 025HR7 Omit Handrim NC

**Projection Handrim**

*NA on plastic coated*

- 025HR8 8 Oblique \$240.00  
*NA with 22" wheels*

**Wheel Lock**

- 025WL1 High Mount Push NC
- 025WL2 High Mount Pull NC
- 025WL5 High Mount Scissor NC  
*NA with grade aids or extension handles*
- 025WL12 Do Not Mount NC
- 025WL10 Omit Wheel Lock NC

**Grade Aid**

- 025GA1 E0974 Grade Aid \$150.00  
*NA with polyurethane, turbo or high-pressure clincher tires*

**Extension Handle** *(NA with Xtender Power Assist)*

- 025EH1 E0961 6" Pair \$75.00

**ACCESSORIES**

**Xtender Power Assist**

- 025XT1 E0986 Xtender Power Assist \$6,595.00  
*Only available with 24" or 26" heat treated wheel rims*

**Xtender Pwr Asst Handrim**

*You must choose one of the following torque options*

- 025XH1 3x Torque Multiplier (Plastic Coated) NC
- 025XH2 1.5x Torque Multiplier (Stainless Steel) NC

**Xtender Pwr Asst Tire**

- 025XW1 E2211 Xtender Pneumatic Tire NC  
E2212

**Xtender Pwr Asst Acc**

- 025XA1 Xtndr Spare NiMH Power Pack \$1,395.00

**Backpack**

- 025A1 Black Utility Backpack \$55.00
- 025A2 Black/Charcoal Kids \$50.00

**Seat Pouch**

- 025A8 Black \$45.00

**Crutch Holder**

*Only available mounted to left side from rear*

- 025A53 E2207 Left Side Mount \$195.00

**Extra Leg Strap**

- 025A55 K0038 Extra Leg Strap Black \$30.00

**Miscellaneous**

- 025A83 Tool Kit \$45.00

**Positioning Belt**

- 025A142 E0978 1.5" Auto Buckle \$70.00
- 025A12 E0978 2" Aircraft Buckle \$110.00
- 025A13 E0978 2" Aircraft Buckle Padded \$120.00

**Spoke Guard**

*NA with Xtender Power Assist*

- 025A23 K0065 Black \$115.00
- 025A24 K0065 Clear \$115.00

**Wheelchair Tray Table**

- 025A31 E0950 Tray Table \$250.00

**Side Guard**

*NA with height adjustable armrests*

- 025SG4 K0108 Plastic Kids \$125.00
- 025SG3 K0108 Plastic Regular \$125.00

**TOTAL CHAIR PRICE \$**

### GP Series Front STF Chart

\*All measurements given are an estimate and have been rounded to the nearest .5"

#### SHORT Frame ( 9 1/4" Front / 8 1/4" Rear)

Fork	Spacer	Fork Hole	3"	4" Poly	5" Poly	6" Poly	6" Pneu	8" Poly	8" Pneu	8x2
				4" Lighted	5" x 1.5	5" Soft Roll				
4 1/4	STD	Bottom	15.50	16.00	16.50	17.00	17.00	x	x	x
4 1/4	3/4	Bottom	16.50	17.00	17.00	17.50	18.00	x	x	x
4 1/4	1 1/2	Bottom	17.00	17.50	18.00	18.50	18.50	x	x	x
5 1/4	STD	Top	x	x	16.50	17.00	17.00	x	x	x
5 1/4	STD	Bottom	16.50	17.00	17.5	18.00	18.00	19.00	20.5	x
5 1/4	3/4	Top	x	x	17.50	18.00	18.00	x	x	x
5 1/4	3/4	Bottom	17.50	18.00	18.50	18.50	19.00	19.50	20.00	x
5 1/4	1 1/2	Top	x	x	18.00	18.50	18.50	x	x	x
5 1/4	1 1/2	Bottom	18.00	18.50	19.00	19.50	19.50	20.50	20.50	x
6	STD	Top	x	x	17.50	18.00	18.00	18.50	19.00	x
6	STD	Bottom	17.50	18.00	18.50	19.00	19.00	19.50	20.00	x
6	3/4	Top	x	x	18.00	18.50	18.50	19.50	19.50	x
6	3/4	Bottom	18.00	18.50	19.00	19.50	19.50	20.50	20.50	x
6	1 1/2	Top	x	x	19.00	19.50	19.50	20.00	20.50	x
6	1 1/2	Bottom	19.00	19.50	20.00	20.50	20.50	21.00	21.50	x
7	STD	Top	x	x	18.50	19.00	19.00	19.50	20.00	x
7	STD	Bottom	18.50	19.00	19.50	20.00	20.00	20.50	21.00	x
7	3/4	Top	x	x	19.00	19.50	19.50	20.50	20.50	x
7	3/4	Bottom	19.00	19.50	20.00	20.50	20.50	21.50	21.50	x
7	1 1/2	Top	x	x	20.00	20.50	20.50	21.00	21.50	x
7	1 1/2	Bottom	20.00	20.50	21.00	21.00	21.50	22.00	22.50	x
8x2	STD	Top	x	x	x	x	x	x	x	19.00
8x2	STD	Bottom	x	x	x	x	x	x	x	20.00
8x2	3/4	Top	x	x	x	x	x	x	x	19.50
8x2	3/4	Bottom	x	x	x	x	x	x	x	20.50
8x2	1 1/2	Top	x	x	x	x	x	x	x	20.50
8x2	1 1/2	Bottom	x	x	x	x	x	x	x	21.00

**STANDARD Frame (10" Front / 9" Rear) ADD 3/4" to each measurement above**

**TALL Frame (11" Front / 10" Rear) ADD 1 3/4" to each measurement above**

\*\*\* Custom Frame Height through Configurator - Front 9 - 13" - Calculate Accordingly (Minimum Front STF = 15 1/2" / Max = 26")

### GP Series Rear STF Chart

\*All measurements given are an estimate and have been rounded to the nearest .5"

SHORT Frame ( 9 1/4" Front / 8 1/4" Rear)			Rear Axle Plate Hole Position						TALL Frame ( 11" Front / 10" Rear)			Rear Axle Plate Hole Position					
Wheel	Rim	Tire	Standard			Inverted			Wheel	Rim	Tire	Standard			Inverted		
			Top	Middle	Bottom	Top	Middle	Bottom				Top	Middle	Bottom	Top	Middle	Bottom
22"	Mag / Spoke	Full Poly	14.5	15	15.5	16	16.5	17	22"	Mag / Spoke	Full Poly	16	16.5	17.5	17.5	18	19
22"	Mag / Spoke	Pneum / Insert	14.5	15	16	16	17	17.5	22"	Mag / Spoke	Pneum / Insert	16.5	17	17.5	18	18.5	19
24"	Mag / Spoke	Full Poly	15	15.5	16.5	16.5	17	17.5	24"	Mag / Spoke	Full Poly	17	17.5	18	18.5	18.5	19.5
24"	Spoke / Spinerger	Clincher	15	16	16.5	16.5	17	18	24"	Spoke / Spinerger	Clincher	17	17.5	18	18.5	19	19.5
24"	Mag / Spoke / Spinerger	Pneum / Insert	15.5	16	17	17	17.5	18	24"	Mag / Spoke / Spinerger	Pneum / Insert	17	17.5	18.5	19	19	20
26"	Spoke / Spinerger	Clincher	16	17	17.5	17.5	18	19	26"	Spoke / Spinerger	Clincher	18	18.5	19	19.5	20	20.5
26"	Spoke / Spinerger	Pneumatic	16.5	17	17.5	18	18.5	19	26"	Spoke / Spinerger	Pneumatic	18	19	19.5	19.5	20	21

STANDARD Frame ( 10" Front / 9" Rear)			Rear Axle Plate Hole Position					
Wheel	Rim	Tire	Standard			Inverted		
			Top	Middle	Bottom	Top	Middle	Bottom
22"	Mag / Spoke	Full Poly	15	15.5	16.5	16.5	17.5	18
22"	Mag / Spoke	Pneum / Insert	15.5	16	17	17	17.5	18
24"	Mag / Spoke	Full Poly	16	16.5	17	17	18	18.5
24"	Spoke / Spinerger	Clincher	16	16.5	17	17.5	18	18.5
24"	Mag / Spoke / Spinerger	Pneum / Insert	16	16.5	17.5	17.5	18	18.5
26"	Spoke / Spinerger	Clincher	17	17.5	18	18.5	19	19.5
26"	Spoke / Spinerger	Pneumatic	17	18	18.5	18.5	19.5	20

Maximum Dump 2 3/4"

\*\*\* Custom Frame Height Through Configurator - REAR 6 1/2 - 13" (Max frame squeeze = 3 1/2") Calculate Accordingly