

QUICKIE V521 POWER (EIPW1, K0011)

July 2004

Account Information	Account No: _____	Drop Ship/Ship To: _____
Date: _____	Name: _____	Address: _____
P.O. #: _____	City: _____	State: _____ Zip: _____
Buyer: _____	Phone: () _____	
Marked For: _____		

QUICKIE V521 POWER FRAMES

QUICKIE V521 POWER			
<input type="checkbox"/>	111PC1	K0011 Quickie V521	\$3,497.00
<input type="checkbox"/>	111PC2	Quickie V521 without Electronics	\$2,199.00
Frame Color			
<input type="checkbox"/>	111C2	Black	NC
<input type="checkbox"/>	111C10	Black Opal	NC
<input type="checkbox"/>	111C1	Blue	NC
<input type="checkbox"/>	111C38	Blue Opal	NC
<input type="checkbox"/>	111C12	Blue Velvet	NC
<input type="checkbox"/>	111C36	Candy Blue	NC
<input type="checkbox"/>	111C35	Candy Red	NC
<input type="checkbox"/>	111C39	Candy Turquoise	NC
<input type="checkbox"/>	111C6	Forest Green	NC
<input type="checkbox"/>	111C4	Midnight Purple	NC
<input type="checkbox"/>	111C37	Orange	NC
<input type="checkbox"/>	111C5	Pearl Pink	NC
<input type="checkbox"/>	111C40	Razzleberry	NC
<input type="checkbox"/>	111C11	Toxic Green	NC
<input type="checkbox"/>	111C8	Yellow	NC
Custom Colors			
<input type="checkbox"/>	111C41	Berry Vibrations	\$250.00
<input type="checkbox"/>	111C20	Graphite	\$250.00
<input type="checkbox"/>	111C42	Stars and Stripes	\$250.00
<input type="checkbox"/>	111C43	Mossy Oak Camouflage	\$250.00

SEATING

Seat Width			
<input type="checkbox"/>	111W5	14" ¹	NC
<input type="checkbox"/>	111W7	16"	NC
<input type="checkbox"/>	111W9	18"	NC
<input type="checkbox"/>	111W11	E2340 20"	\$325.00
Seat Depth			
<input type="checkbox"/>	111SD7	14"	NC
<input type="checkbox"/>	111SD9	16"	NC
<input type="checkbox"/>	111SD11	18"	NC
Seat Height			
<input type="checkbox"/>	111HT11	18-1/2"	NC
<input type="checkbox"/>	111HT15	20-1/4"	NC
Seat Sling			
<input type="checkbox"/>	111SL4	Seat Sling	NC
<input type="checkbox"/>	111SL6	Omit	NC
Cushion			
<input type="checkbox"/>	111CU1	E0963 2" Cushion	\$95.00
Back Post Type			
<input type="checkbox"/>	111BP1	Straight Back Post	NC
<input type="checkbox"/>	111BP3	8° Bend Back Post	\$25.00
Back Type			
<input type="checkbox"/>		Fixed Back	NC
<input type="checkbox"/>		Angle Adjustable Back	\$165.00

Back Upholstery

<input type="checkbox"/>	111BU1	Back Upholstery Standard	NC
<input type="checkbox"/>	111BU14	Omit Back Upholstery	-\$15.00
<input type="checkbox"/>	111BP8	Omit Back Post Upholstery	-\$35.00

Back Height

<input type="checkbox"/>	111BP35	15"	NC
<input type="checkbox"/>	111BP36	16"	NC
<input type="checkbox"/>	111BP37	17"	NC
<input type="checkbox"/>	111BP38	18"	NC
<input type="checkbox"/>	111BP39	19"	NC
<input type="checkbox"/>	111BP40	20"	NC

ARMREST

Armrest			
<input type="checkbox"/>	111AA25	Dual Post Flip Back	NC
<input type="checkbox"/>	111AA42	E0973 Dual Post Flip Back Height Adjustable	\$205.00
<input type="checkbox"/>	111AA8	E0973 Single Post Height Adjustable	\$195.00
<input type="checkbox"/>	111AA14	Omit Armrest	NC

Armpad

<input type="checkbox"/>	111AM1	Desk Length (10")	\$15.00
<input type="checkbox"/>	111AM2	Full Length (14")	\$15.00
<input type="checkbox"/>	111AM7	Desk Length Waterfall (10")	NC
<input type="checkbox"/>	111AM8	Full Length Waterfall (14")	NC

FOOTREST

Hanger			
<input type="checkbox"/>	111H2	70° Swing Away	NC
<input type="checkbox"/>	111H37	90° Swing Away ²	NC
<input type="checkbox"/>	111H14	E0990 Elevating ^{3,4}	\$220.00
<input type="checkbox"/>	111H32	Omit Hanger	-\$15.00

Extension Tube ⁵

<input type="checkbox"/>	111E1	Short ⁴	NC
<input type="checkbox"/>	111E3	Medium	NC
<input type="checkbox"/>	111E18	Long ⁴	NC

Footplate

<input type="checkbox"/>	111F1	Composite w/Heel Loop	NC
<input type="checkbox"/>	111F3	K0040 Adult Angle Adjustable w/Heel Loop	\$105.00
<input type="checkbox"/>	111F9	K0040 Kids Angle Adj. W/Leg Strap ⁶	\$105.00
<input type="checkbox"/>	111F47	K0040 Multi-position w/Heel Loop	\$150.00
<input type="checkbox"/>	111F5	Aluminum Footplate ⁷	\$30.00
<input type="checkbox"/>	111F6	Omit Footplate	-\$15.00

Footrest Accessories

<input type="checkbox"/>	111A110	E0952 Toe Loop	\$20.00
<input type="checkbox"/>		K0038 Leg Strap (Black)	\$30.00

FOOTNOTES:

- 1 Only available with U1 batteries.
- 2 Only available with kids angle adjustable footplates.
- 3 Requires aluminum or angle adjustable footplate.
- 4 Not available with multi-position footplate.
- 5 Not available with 90° hangers.
- 6 Only available with 90° hangers.
- 7 Only available with elevating legrest and does not come with heel loops.

Note: The HCPCS codes provided are based on code verification by the SADMERC or our interpretation of the code definitions.

Proper use of HCPCS codes does not ensure coverage or payment. For coverage information, verify the appropriate payer's coverage policy.

CUSTOMER SERVICE

Phone: (800) 236-4215
Fax: (602) 230-9497

ELECTRONICS

Joystick			
<input type="checkbox"/>	111J20	Integral.....	NC
Mount Type			
<input type="checkbox"/>	111JM1	Fixed Mount.....	NC
<input type="checkbox"/>	111JM2	E1028 Swing Away Mount.....	\$235.00
Hand Mount			
<input type="checkbox"/>	111HM1	Right.....	NC
<input type="checkbox"/>	111HM2	Left.....	NC
Height Adjustable Mount			
<input type="checkbox"/>		Fixed Mount.....	NC
<input type="checkbox"/>		Adjustable Mount.....	\$95.00
Joystick Handle			
<input type="checkbox"/>	111JH3	Stick.....	NC
<input type="checkbox"/>	111JH4	E2323 T-Handle.....	\$75.00
<input type="checkbox"/>	111JH1	E2323 Ball.....	\$75.00
<input type="checkbox"/>	111JH2	E2323 Foam Ball.....	\$75.00
Battery Charger			
<input type="checkbox"/>		5 A ¹	NC
<input type="checkbox"/>		8 A Dual Mode ²	\$85.00
<input type="checkbox"/>	111BJ3	Omit Battery Charger.....	-\$40.00
Battery Box			
<input type="checkbox"/>	111BI4	1 Box for U1 Batteries.....	NC
<input type="checkbox"/>	111BI5	2 Box for 22NF Batteries ^{3,4}	\$75.00
DRIVE WHEEL			
Drive Wheels			
<input type="checkbox"/>	111DW10	12" Pneumatic Tire.....	NC
<input type="checkbox"/>	111DW7	K0093 12" Pneumatic w/Insert.....	\$220.00
<input type="checkbox"/>	111DW12	12" Solid.....	\$110.00
Wheel Lock			
<input type="checkbox"/>	111WL8	Push to Lock.....	\$85.00
<input type="checkbox"/>	111WL12	Do Not Mount.....	NC
Wheel Lock Extension Handle			
<input type="checkbox"/>	111EH1	Extension Handle 6".....	\$50.00
CASTERS			
Caster Wheel			
<input type="checkbox"/>	111CW21	8" Solid.....	NC
<input type="checkbox"/>	111CW23	8" Pneumatic.....	\$85.00
<input type="checkbox"/>	111CW22	K0097 8" Pneumatic w/Insert.....	\$125.00
ACCESSORIES			
Batteries ^{3,5}			
<input type="checkbox"/>	111A46	E2365 Battery Gel Cel U1 (pair).....	\$230.00
<input type="checkbox"/>	111A47	E2361 Battery Gel Cel 22NF (pair) ^{4,5,6}	\$310.00
Backpack			
<input type="checkbox"/>	111A1	Black Utility Backpack.....	\$55.00
<input type="checkbox"/>	111A2	Black/Charcoal Kids.....	\$50.00
Holders			
<input type="checkbox"/>	K0104	Oxygen Tank Holder.....	\$150.00

Seat Pouch

<input type="checkbox"/>	111A8	Black.....	\$45.00
Miscellaneous			
<input type="checkbox"/>	111A82	Caddy.....	\$45.00
Positioning Belt			
<input type="checkbox"/>	111A12	E0978 Belt 2" Aircraft Buckle.....	\$85.00
<input type="checkbox"/>	111A13	E0978 Belt 2" Aircraft Buckle Padded.....	\$115.00
<input type="checkbox"/>		E0978 Belt 1.5" Auto Buckle.....	\$50.00
Touch Up Paint			
<input type="checkbox"/>	111A27	Frame Touch Up Paint.....	\$25.00
<input type="checkbox"/>	111A28	X-Tube.....	\$25.00
		Color: _____	
Wheelchair Tray Table ⁷			
<input type="checkbox"/>	111A36	E0950 Medium (14"-16").....	\$225.00
<input type="checkbox"/>	111A37	E0950 Large (17"-20").....	\$225.00
Programmer			
<input type="checkbox"/>	111PR1	Qtronix Programmer Pad.....	\$405.00

TOTAL CHAIR PRICE \$ _____

Cushion _____

Back _____

ADDITIONAL QUICKIE V521 SPECIFICATIONS

Specifications subject to change without notice.

Weight: 83 lbs. (Without batteries & Front Riggings)

Overall Width: 19.5" (14" Frame Width)
 21.5" (16" Frame Width)
 23.5" (18" Frame Width)
 25.5" (20" Frame Width)

Overall Length: 30" (14" Seat Depth & Without Front Riggings)
 32" (16" Seat Depth & Without Front Riggings)
 34" (18" Seat Depth & Without Front Riggings)

Top Speed: 4.5 mph

FOOTNOTES:

- 1 _____
- 2 Recommended for gel cell and sealed lead acid batteries only, Use 8 A dual mode charger for wet cell batteries.
- 3 Recommended for wet cell batteries.
- 4 Recommended for wet cell batteries.
- 5 Separate battery boxes accommodates either U1 or 22NF battery sizes.
- 6 Not available with 14" seat width.
- 7 22NF batteries recommended for best range. Must use 2 battery boxes for 22NF batteries. Requires standard armpads to mount.

CUSTOMER SERVICE

Phone: (800) 236-4215
 Fax: (602) 230-9497