



Quickie XTR

EQS, Effective August 7, 2006

ACCOUNT INFORMATION

Account No: _____
 Date: _____
 P.O. #: _____
 Buyer: _____
 Marked For: _____

Drop Ship/Ship To: _____
 Name: _____
 Address: _____
 City: _____ State: _____
 Zip Code: _____ Phone Number: _____

QUICKIE XTR

K0005 Quickie XTR \$1,999.00

FRAMES

Shock Size

089SK1 Light Shock (Up to 160 lbs) ¹ NC
 089SK2 Heavy Shock (Up to 250 lbs) ¹ NC

Frame Color

089C2 Black NC
 089C10 Black Opal NC
 089C1 Blue NC
 089C38 Blue Opal NC
 089C12 Blue Velvet NC
 089C36 Candy Blue NC
 089C35 Candy Red NC
 089C39 Candy Turquoise NC
 089C6 Forest Green NC
 089C4 Midnight Purple NC
 089C37 Orange NC
 089C5 Pearl Pink NC
 089C40 Razzleberry NC
 089C11 Toxic Green NC
 089C8 Yellow NC

Custom Colors

089C41 Berry Vibrations \$210.00
 089C20 Graphite \$210.00
 089C15 Purple Chase \$210.00
 089C42 Stars and Stripes \$210.00

Frame Width

089W5 14" NC
 089W6 15" NC
 089W7 16" NC
 089W8 17" NC
 089W9 18" NC
 089W10 19" \$150.00
 089W11 E2201 20" \$300.00

SEATING

Sling Depth

089SD7 14" NC
 089SD8 15" NC
 089SD9 16" NC
 089SD10 17" NC
 089SD11 18" NC

Cushion

089CU1 2" Cushion \$95.00
 089CU2 Omit Cushion NC
 089CU6 Omit Seat Sling -\$15.00

Back Type

089BT1 Folding Backrest NC

Back Height

089B78 11" - 14" NC
 089B83 12" - 16" NC
 089B94 16" - 20" NC

Back Upholstery

089BU1 Standard NC
 089BU8 Adjustable ² \$150.00
 089BU16 Reinforced ³ \$30.00
 089BU14 Omit Upholstery -\$15.00
 089BT14 Omit Backrest & Upholstery -\$35.00

Back Rest Options

089BT33 Folding Backrest with Integrated Push Handles NC
 089A85 Bolt-On Push Handles \$75.00

FOOTREST

Footrest

089F57 Rigid Footrest Standard NC
 089F50 Platform Footrest NC

Adaptive Platform Spacer ⁴

089AD1 2" \$95.00
 089AD2 4" \$95.00
 089AD3 6" \$95.00

Footrest Accessories

089LS6 K0038 Leg Strap NC
 089IG2 Neoprene Impact Guard \$85.00

ARMREST

Armrest

089AA15 Padded Swing-Away Armrest \$175.00

CASTERS

Caster

089CW1 3" Micro NC
 089CW30 3" Light Up \$50.00
 089CW3 4" Low Profile Polyurethane Tire ⁵ NC
 089CW5 5" Low Profile Polyurethane Tire NC

Suspension Fork

089CF33 E1015 Frog Leg Caster Fork MC5 ⁶ \$330.00

AXLE AND WHEEL

Camber Tube

089CB1 0° Camber NC
 089CB3 4° Camber NC
 089CB5 8° Camber NC

FOOTNOTES:

- ¹ Recommended body weight option.
- ² Not available with fabric side guards.
- ³ Not available with adjustable upholstery.
- ⁴ Only available with platform footrest.
- ⁵ 4" Low Profile Polyurethane only available with Frog Legs MC5 Suspension Forks.
- ⁶ Only available with 4" Casters.

PROPER USE OF HCPCS CODES DOES NOT ENSURE COVERAGE OR PAYMENT. FOR COVERAGE INFORMATION, VERIFY THE APPROPRIATE PAYER'S COVERAGE POLICY.

Extra Camber		
<input type="checkbox"/>	089XB1	0° Camber \$35.00
<input type="checkbox"/>	089XB3	4° Camber \$35.00
<input type="checkbox"/>	089XB5	8° Camber \$35.00

Axles		
<input type="checkbox"/>	089AX4	Quick Release Steel NC
<input type="checkbox"/>	089AX5	Quick Release Titanium ² \$135.00
<input type="checkbox"/>	089AX11	Quad Quick Release \$90.00
<input type="checkbox"/>	089AX16	Quad Quick Release Titanium ² \$225.00

Rear Anti Tip		
<input type="checkbox"/>	089AT1	E0971 Rear Anti-Tip \$135.00
<input type="checkbox"/>	089AT3	Omit Rear Anti-Tip NC

Rear Wheel		
<input type="checkbox"/>	089RW1	24" Mag (540) ³ NC
<input type="checkbox"/>	089RW2	24" Spoke (540) ³ NC
<input type="checkbox"/>	089RW11	26" Spoke (590) ^{4,5} NC
<input type="checkbox"/>	089RW19	24" QP Heat-Treated Basketball (559) \$250.00
<input type="checkbox"/>	089RW23	24" QP Heat-Treated Tennis 32H (559) \$250.00
<input type="checkbox"/>	089RW32	24" Mountain (540) Knobbie Tires Only \$250.00
<input type="checkbox"/>		24" Spinergy (540)(Clincher Only) ^{8,9} \$695.00
<input type="checkbox"/>	089RR26	Black Spokes
<input type="checkbox"/>	089RR15	Yellow Spokes
		26" Spinergy (Clincher Only) ^{8,9} \$695.00
<input type="checkbox"/>	089RR26	Black Spokes
<input type="checkbox"/>	089RR15	Yellow Spokes
<input type="checkbox"/>	089RR25	Omit Rear Wheel and Axle -\$65.00
		(Intended size for set up)

Rear Wheel Tire		
<input type="checkbox"/>	089RT1	Pneumatic ¹⁰ NC
<input type="checkbox"/>	089RT2	K0064 Pneumatic Airless Insert ^{10,11} \$70.00
<input type="checkbox"/>	089RT6	High Pressure Clincher ¹⁰ \$95.00
<input type="checkbox"/>	089RT8	Kevlar ^{10,11} \$95.00
<input type="checkbox"/>	089RT5	Turbo Tread ¹² NC

Handrim		
<input type="checkbox"/>	089HR1	Aluminum Anodized NC
<input type="checkbox"/>	089HR4	Plastic Coated \$85.00
<input type="checkbox"/>	089HR7	Omit Handrim -\$20.00

Projections (NA of Plastic Coated)		
<input type="checkbox"/>	089HR8	8 Oblique Handrim ¹³ \$195.00
<input type="checkbox"/>	089HR9	8 Vertical Handrim ¹³ \$195.00
<input type="checkbox"/>	089HR10	12 Oblique Handrim ^{12,13} \$195.00
<input type="checkbox"/>	089HR11	12 Vertical Handrim ^{12,13} \$195.00

Wheel Lock		
<input type="checkbox"/>	089WL1	High Mount Push NC
<input type="checkbox"/>	089WL10	Omit Wheel Lock NC

ACCESSORIES

Back Pack		
<input type="checkbox"/>	089A1	Black Utility Backpack \$55.00
<input type="checkbox"/>	089A2	Black/Charcoal Kids \$50.00

Seat Pouch		
<input type="checkbox"/>	089A8	Black \$45.00

Extra Leg Strap		
<input type="checkbox"/>	K0038	Extra Leg Strap \$30.00

Miscellaneous		
<input type="checkbox"/>	089A86	Basketball Roller \$45.00
<input type="checkbox"/>	089A83	Tool Kit \$25.00

Positioning Belt			
<input type="checkbox"/>	E0978	1.5" Auto Buckle	\$60.00
<input type="checkbox"/>	089A12	E0978 2" Aircraft Buckle	\$85.00
<input type="checkbox"/>	089A13	E0978 2" Aircraft Buckle Padded	\$115.00

Spoke Guard			
<input type="checkbox"/>	089A23	K0065 Black ²	\$90.00
<input type="checkbox"/>	089A24	K0065 Clear ²	\$90.00
<input type="checkbox"/>	089A25	K0065 Sport Hi Profile ²	\$120.00
<input type="checkbox"/>	089A26	K0065 Sport Low Profile ²	\$120.00

Touch Up Paint ¹⁴		
<input type="checkbox"/>	089A27	Frame \$25.00

Side Guard		
<input type="checkbox"/>	089SG1	Fabric Regular \$45.00
<input type="checkbox"/>	089SG4	Plastic Kids \$120.00
<input type="checkbox"/>	089SG3	Plastic Regular \$120.00

TOTAL CHAIR PRICE

\$ _____

JAY SEATING OPTIONS*

Select Cushion	Select Back
<input type="checkbox"/> Basic (E2601/E2602)	<input type="checkbox"/> Basic (E2611)
<input type="checkbox"/> Soft Combi P (E2605, E2606)	<input type="checkbox"/> J2 (E2613)
<input type="checkbox"/> Xtreme (E2603)	<input type="checkbox"/> J2 Deep Contour (E2615)
<input type="checkbox"/> Xtreme S (E2603)	<input type="checkbox"/> Xtreme
<input type="checkbox"/> Triad (E2603)	<input type="checkbox"/> J2 Plus (E2613/E2614)
<input type="checkbox"/> Triad P (E2603)	<input type="checkbox"/> Care (E2615)
<input type="checkbox"/> Care (K0669)	
<input type="checkbox"/> Duo (E2603)	
<input type="checkbox"/> Easy (E2607)	
<input type="checkbox"/> Jay 2 (K0108)	
<input type="checkbox"/> J2 Recline (K0108)	
<input type="checkbox"/> J2 Deep Contour (K0108)	

Indicate Backrest Width

Indicate Exact Cushion Dimensions (w" x d") _____

If a modification is necessary to achieve your cushion or backrest dimensions, do you want customer service to notify you prior to placing this order?

Please Choose

- Send Quote for any Modification
- Send Quote only when a modification has a charge
 - Fax me # _____
 - Email me _____

* Please see Jay Price List for latest pricing information. Also, note that Jay cushions and backs may ship from different distribution centers.

FOOTNOTES:

- ¹ Only available on heat-treated wheels.
- ² Not available with Spinergy wheels.
- ³ Only available with standard pneumatic tire and Kevlar.
- ⁴ Not available with heat-treated wheels.
- ⁵ Only available with standard pneumatic tire or high pressure clincher tire.
- ⁶ Only available with Turbo Tread tires.
- ⁷ Not considered everyday or indoor tire.
- ⁸ Not available with spoke guards.
- ⁹ Only available with yellow or black spokes.
- ¹⁰ Not available with performance wheels.
- ¹¹ Only available on 24" wheels.
- ¹² Heat-treated wheels only.
- ¹³ Not available with Mountain or Spinergy wheels.
- ¹⁴ Not available with Custom Colors.

HIGH PRESSURE TIRES

	0 Degree Camber Tube		4 Degree Camber Tube		8° Degree Camber Tube	
	Front Seat Height	Rear Seat Height	Front Seat Height	Rear Seat Height	Front Seat Height	Rear Seat Height
14" Depth	18.75	17.00	18.75	16.50	18.25	16.00
15" Depth	19.00	17.00	19.00	16.50	18.25	16.00
16" Depth	19.00	17.00	19.00	16.50	18.50	16.00
17" Depth	19.25	17.00	19.25	16.50	18.50	16.00
18" Depth	19.25	17.00	19.25	16.50	18.75	16.00

** rounded to nearest 1/4" * with 8° camber, caster is mounted in upper hole set