



Quickie Shark

EISP5, Effective July 2009

ACCOUNT INFORMATION

Account No: _____
 Date: _____
 P.O. #: _____
 Buyer: _____
 Marked For: _____

Drop Ship/Ship To: _____
 Name: _____
 Address: _____
 City: _____ State: _____
 Zip Code: _____ Phone Number: _____

Quickie Shark

■ 157BX8 Quickie Shark \$3,349.00

FRAMES

Frame Color

- 157C91 Amethyst Fantasy \$95.00
- 157C55 Black NC
- 157C1 Blue NC
- 157C13 Candy Teal NC
- 157C75 Flip Flop Blue / Green \$95.00
- 157C61 Forrest Green NC
- 157C85 Graphmetallic NC
- 157C92 Illusion Lapis \$95.00
- 157C93 Mystery \$95.00
- 157C84 Platinum metallic NC
- 157C54 Red Ruby NC
- 157C90 Silver NC
- 157C16 Silver Texture NC
- 157C53 Sun Yellow NC
- 157C94 Tutti Frutti Orange NC
- 157C95 White RAL NC
- 157C34 Omit Color Frame NC

Brilliant Colors:

- 157C73 Brilliant Lime NC
- 157C72 Brilliant Marine NC
- 157C52 Brilliant Purple NC
- 157C74 Brilliant Red NC
- 157C76 Brilliant Silver NC

Matte Colors:

- 157C68 Matte Lime NC
- 157C96 Matte Candy Teal NC
- 157C67 Matte Marine NC
- 157C65 Matte Purple NC
- 157C63 Matte Red NC
- 157C62 Matte Silver NC

Frame Length

- 157FL3 Frame Length Short NC

Frame Length Notes:


Overall with:
 20" Wheels 77.56" (197 cm)
 26" Wheels 79.53" (202 cm)
 without Rear Wheels 55.91" (142 cm)
 without steering fork & wheels 41.34" (105 cm)
 Transfer Seat Height 8.66" (22 cm)

Position Fork Mount

- 157PZ01 15" Setting (88cm) NC
 - 157PZ02 15.75" Setting (40 cm) NC
 - 157PZ03 16.5" Setting (42 cm) NC
- w/adjustable ball bearing

Position Fork Mount Diagram:

High Mount Fork Preset
 @ 15.75" (40 cm)



Fork Style

- 157FK03 Fixed Fork (standard frame) NC
- 157FK01 Detachable Fork \$250.00
- 157FK02 X-Link Fork (removable) \$450.00

X-Link Fork Option

Easier access-
 Raises seat 3.94" (10 cm) higher for transfers
 Moves the crank 5.91" (15 cm) forward for transfers

CONTROL OPTIONS

Hand Crank

- 157CK1 Crank Standard V Crank NC

Handle Type

- 157HY4 Horizontal Handle^{1,3} NC
- 157HY6 Vertical Handle Small⁴ NC
- 157HY7 Vertical Handle Std⁴ NC

Hand Mount

- 157HM5 On Frame Right Side² NC
- 157HM6 Handle Right Side NC

Control Notes:

Standard V Crank
 6.89" (17.5 cm) long
 19.29" (49 cm) wide (center to center of Handle)
Handle
 Vertical Small and Std @ 15 Degree Diagonal
Brake
 Parking Brake Std
 Brake lever can be locked when brake is on

FOOTNOTES:

- ¹ Only available with 7 Speed Option.
- ² Not available with Detachable or X-Link Fork Style.
- ³ Horizontal Handle with 7 Speed only available with Steering Fork Right Side Hand Mount
- ⁴ NA with frame mounted right side.

PROPER USE OF HCPCS CODES DOES NOT ENSURE COVERAGE OR PAYMENT. FOR COVERAGE INFORMATION, VERIFY THE APPROPRIATE PAYER'S COVERAGE POLICY.

Speed		
<input type="checkbox"/>	157SA4	3 X 9 Gear (27 speed) NC

Speed Notes:
 7 Speed with integral back-pedal brake
 27 Speed with rim brake on left handle

SEATING


Seat Sling		
<input type="checkbox"/>	157SL1	Seat Sling Standard NC

Seat Width		
<input type="checkbox"/>	157S6	15" Seat Width (38 cm) NC
<input type="checkbox"/>	157S21	16.5" Seat Width (42 cm) NC
<input type="checkbox"/>	157S9	18" Seat Width (46cm) NC

BACKREST


Back Type		
<input type="checkbox"/>	157BT10	Standard Backrest NC
		<i>Adjustable height & angle adjustable</i>
<input type="checkbox"/>	157BT1	Folding Backrest \$140.00
<input type="checkbox"/>		Medium height and angle adjustable with quick release NC

Backrest Notes:
 Standard Backrest (shown on right)
 22.83" (58 cm) long
 Height and Angle Adjustable
 (30 to 90 degree)
Fold Down Backrest
 22.83" (58 cm) long
 Uses Quick Release Pin and Retainer Ring for ease of adjustment




Backrest Position		
<input type="checkbox"/>	157BQ07	Adj (10"-21.5") Preset @ 12" NC

Backrest Position Diagram:
 Preset @ 12" (30.48 cm)



Neck Support		
<input type="checkbox"/>	157NT01	Neck Support \$125.00

Neck Support Diagram:
 Height and Depth Adjustable



FOOTREST

Footrest		
<input type="checkbox"/>	157F236	Adj Short (16"-31.5") (40-80 cm) NC
<input type="checkbox"/>	157F237	Adj Long (27.5"-37.5") (70-95 cm) NC

Footrest Notes:
 Adjustable Short (16"-31.5")
 Adjustable Long (27.5"-37.5")

WHEELS

Drive Wheel		
<input type="checkbox"/>	157DW33	26 X 1.35" Marathon NC
<input type="checkbox"/>	157DW34	26 X 1" Hi Perf \$75.00
<input type="checkbox"/>	157DW35	Omit Drive Wheel -\$50.00

Rear Wheel Size		
<input type="checkbox"/>	157WS4	20" Wheel NC
<input type="checkbox"/>	157WS18	26" Wheel Dt Rader \$275.00
<input type="checkbox"/>		Omit Rear Wheels -\$55.00

Rear Wheel Tire		
<input type="checkbox"/>	157RT19	1.35" Marathon Tire NC
<input type="checkbox"/>	157RT20	1" Hi Performance Reifen \$90.00

Axle		
<input type="checkbox"/>	157AX17	Threaded Axle NC
<input type="checkbox"/>	157AX01	Quick Release Axle NC

ACCESSORIES

Accessories Sports		
<input type="checkbox"/>	157AZ01	Front Fender (26") \$50.00
<input type="checkbox"/>	157AZ02	Rear Fender (20") \$75.00
<input type="checkbox"/>		Rear Fender (26") \$75.00
<input type="checkbox"/>	157AZ03	Front Headlight \$50.00
<input type="checkbox"/>	157AZ04	Rear Light \$50.00
<input type="checkbox"/>	157AZ05	Cushion W/Composite Foam \$125.00
<input type="checkbox"/>	157AZ06	Waist Belt \$125.00
<input type="checkbox"/>	157AZ07	Signal Flag NC
<input type="checkbox"/>	157AZ08	Rear Bumper 20" Rear Wheels \$150.00
<input type="checkbox"/>	157AZ09	Rear Bumper 26" Rear Wheels \$150.00

HANDCYCLE TOTAL _____

Handling Charge ²

<input type="checkbox"/>	Standard Handling	NC
<input type="checkbox"/>	Expedited Handling	\$895.00

TOTAL _____

FOOTNOTES:

- ¹ Must select Axle type
- ² Standard Handling is a combination of sea & ground freight
 Expedited Handling is a combination of air & ground freight
 Please contact Customer Service for manufacturing lead times.