



V-121



Now standard
with the VSI
joystick controller

Photo may vary from actual product.





V-121

The Quickie V-121 offers convenience and transportability. The V-121 has excellent performance on both indoor and hard outdoor surfaces. The quick-release, removable battery cage makes disassembly and transportation easy. The V-121 is a lightweight, easy-to-fold and easy-to-transport power chair.



Now standard with VSI joystick controller



V-121 Specifications:

Frame Style	Folding Frame
Overall	
Weight	64 lbs without batteries
Overall Width	Seat width + 5.5"
Overall Length (16" seat depth)	31"
Seating	
Seat Widths	14", 16", 18", 20"
Seat Depths	14", 16" & 18"
Seat Height	18.5" & 20"
Back Height	15", 16", 17", 18", 19", 20"
Back Angle	Fixed or adjustable (-5° to +5°) (option)
Positioning Belt	Auto style buckle
Seat & Back Material	Nylon
Footrests	
Hangers	Swing-away & Elevating (option)
Footplates	Flip-up with heel loops & Angle Adjustable (option)
Adjustability	Length adjustable
Armrests	
Style	Flip-back, Height adjustable & Removable
Armpad	Desk & Full Length
Performance	
Maximum user weight	250 lbs
Minimum turning radius*	33.5"
Maximum speed	4.5 mph
Range	20 miles
* Measured using 16"x16" with standard footrests	
Wheels & Casters	
Drive Wheel	Split rim
Drive Tire	12"x2.25" Solid Insert & Pneumatic
Caster Size	8"x2"
Caster Type	Solid & Pneumatic (option)
Anti-tips	Standard
Electronics	
Joystick/Controller	P&G VSI Proportional Integral
Output	50 Amps + 20Amp Boost
Programmable	Speed - Forward, Reverse, Turning Acceleration - Forward, Reverse, Turning (Sensitivity) Joystick Throw, Steer Correct
Drive Profiles	up to 5 Independent Drive Profiles
Battery Charger	24V, 5 Amp off-board charger (wet cell, SLA & gel cell)
Batteries (Option)	U1 or 22NF capacity
Additional Options	
	Oxygen Tank Holder
	Swing-away Joystick Mount
	Programming Pad
	Height Adjustable Joystick Mount
	Backpack
	Wheelchair Tray Table

Quickie V-121 Highlights:

- Excellent performance on indoor and hard outdoor surfaces
- Standard flip-back armrests, optional height-adjustable, flip-back armrest
- Standard swing-away footrests, optional elevating legrest
- 12" drive wheels
- Choice of 14", 16", 18" or 20" seat widths
- Direct-drive motors and gear boxes
- Charger jack/programming port located on joystick for convenient connections
- Optional 22NF batteries for greater range